

Abaxis ASTM Protocol Sample Records

The following sample records are examples of the Abaxis implementation of the ASTM 1394-97 standard within the ASTM 1381-02 low level communications protocol standard. The examples address run time rotor results messages, comment messages, query of rotor result messages, and exception comment messages.

***** Example (Abaxis to LIS) sign on comment message *****

***** piccolo xpress version # 2.1.9, serial # 0000V0000D, at 2006 04/01 06:24:21 *****

```
<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401062421<CR><ETX>A1<CR><LF>
<STX>2C|1|I|^Power Up Initialization Sign On|G<CR><ETX>94<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

***** Example (Abaxis to LIS) sign on comment message *****

***** VetScan VS2 version # 2.1.9, serial # 0000V0000D, at 2006 04/01 06:24:21 *****

```
<STX>1H|¥^&|||ABAXIS, INC.^VetScan VS2^2.1.9^0000V0000D|||||P|E 1394-97|20060401062421<CR><ETX>A1<CR><LF>
<STX>2C|1|I|^Power Up Initialization Sign On|G<CR><ETX>94<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

***** Example (Abaxis to LIS) rotor abort error results message *****

***** patient ID 12345678901234, MALE, ROTOR UNKNOWN, operator ID 12345678901234, doctor ID 12345678901234 *****

```
<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401062828<CR><ETX>AC<CR><LF>
<STX>2P|1|12345678901234|44444|33333|||19530225|M|||789-012-3456|88888^53 Yrs.^Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1||||20060401062445|||12345678901234||||12345678901234|||||F<CR><ETX>93<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

***** Example (Abaxis to LIS) rotor abort error results message *****

***** patient ID 98765432109876, FEMALE, ROTOR UNKNOWN, operator ID 98765432109876, doctor ID 98765432109876 *****

```
<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401063227<CR><ETX>A6<CR><LF>
<STX>2P|1|98765432109876|44444|33333|||19530225|F|||789-012-3456|88888^53 Yrs.^Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1||||20060401062852|||98765432109876||||98765432109876|||||F<CR><ETX>BD<CR><LF>
```

```
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

***** Example (Abaxis to LIS) rotor abort error results message *****

***** patient ID LEVEL II, CONTROL, ROTOR UNKNOWN, operator ID 87654321098765, doctor ID 87654321098765 *****

```
<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401063647<CR><ETX>AC<CR><LF>
<STX>2P|1|ID LEVEL II|||||U|||||^|^Control|||||<CR><ETX>21<CR><LF>
<STX>30|1|||||20060401063314|||||87654321098765|||||87654321098765|||||F<CR><ETX>AF<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

***** Example (Abaxis to LIS) rotor abort error results message *****

***** patient ID 76543210987654, ROTOR UNKNOWN, operator ID 76543210987654, doctor ID 76543210987654 *****

```
<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401064223<CR><ETX>A3<CR><LF>
<STX>2P|1|76543210987654|||||U|||||^|^|||||<CR><ETX>21<CR><LF>
<STX>30|1|||||20060401063833|||||76543210987654|||||76543210987654|||||F<CR><ETX>AD<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

*** Example (LIS to Abaxis) Query message ***

*** patient ID all, all rotors ran have matching query results dumped below ***

<STX>1H|¥&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|all<CR><ETX>F5<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥&||ABAXIS, INC. ^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401065529<CR><ETX>AD<CR><LF>
<STX>2P|1|12345678901234|44444|33333|||19530225|M|||789-012-3456||88888^53 Yrs. ^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1||||20060401062445|||12345678901234|||||12345678901234|||||||Q<CR><ETX>9E<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>

<STX>1H|¥&||ABAXIS, INC. ^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401065532<CR><ETX>A7<CR><LF>
<STX>2P|1|98765432109876|44444|33333|||19530225|F|||789-012-3456||88888^53 Yrs. ^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1||||20060401062852|||98765432109876|||||98765432109876|||||||Q<CR><ETX>C8<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>

<STX>1H|¥&||ABAXIS, INC. ^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401065534<CR><ETX>A9<CR><LF>
<STX>2P|1|ID LEVEL II|||||U|||||^|^Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1||||20060401063314|||87654321098765|||||87654321098765|||||||Q<CR><ETX>BA<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>

```
<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401065536<CR><ETX>AB<CR><LF>
<STX>2P|1|76543210987654|||||U|||||^~|||||<CR><ETX>21<CR><LF>
<STX>30|1|||||20060401063833|||||76543210987654|||||76543210987654|||||Q<CR><ETX>B8<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

***** Example (LIS to Abaxis) Query message *****

***** patient ID 12345678901234, only rotor with matching patient ID query results sent below *****

```
<STX>1H|¥^&|||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|12345678901234<CR><ETX>93<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>
```

```
<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401065609<CR><ETX>AC<CR><LF>
<STX>2P|1|12345678901234|44444|33333|||19530225|M|||789-012-3456|88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||||20060401062445|||||12345678901234|||||12345678901234|||||Q<CR><ETX>9E<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

***** Example (LIS to Abaxis) Query message *****

***** patient ID 87654321098765, no rotor with matching patient ID exists so no query results sent below *****

```
<STX>1H|¥^&|||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|87654321098765<CR><ETX>A3<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>
```

```
<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||P|E 1394-97|20060401065710<CR><ETX>A5<CR><LF>
<STX>2P|1|87654321098765<CR><ETX>A2<CR><LF>
<STX>30|1|||||Z<CR><ETX>39<CR><LF>
<STX>4C|1|I|^~87654321098765:::No Results for this Query|G<CR><ETX>CD<CR><LF>
<STX>5L|1|I<CR><ETX>03<CR><LF>
```

*** patient ID level ii, rotor with matching patient ID level ii exists so query results sent below ***

```
<STX>1H|¥&|||ABAXIS, INC. ^pico xpress^2.1.9^0000V0000|||||P|E 1394-97|20060401065826<CR><ETX>AD<CR><LF>
<STX>2P|1|ID LEVEL II|||||U|||||^~|Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1|||20060401063314|||87654321098765|||87654321098765|||||||Q<CR><ETX>BA<CR><LF>
<STX>4C|1|I|^46FF FFFF^Run Count¥^7DFF FFFF^Abort Count|I<CR><ETX>CF<CR><LF>
<STX>5C|2|I|^28FF FFFF^Print Count¥^5D19 6FFF^Flash Count|I<CR><ETX>9D<CR><LF>
<STX>6C|3|I|^0000^DAC Trim Flags¥^0000^Offset Error Flags¥^0000 0000^System Flags|I<CR><ETX>0B<CR><LF>
<STX>7C|4|I|^0000 0000^Bead Check 1 Flags¥^0000 0000^Bead Check 2 Flags|I<CR><ETX>E5<CR><LF>
<STX>0C|5|I|^0000 0000^Empty Cuvette Flags¥^0000 0000^Distribution Check Flags|I<CR><ETX>75<CR><LF>
<STX>1C|6|I|^0000^Quality Control Flags¥^0000^Offset SD Error Flags¥^0000^Wavelength CV Flags|I<CR><ETX>8F<CR><LF>
<STX>2C|7|I|^4055^BARCODE ERROR|I<CR><ETX>62<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

*** piccolo xpress software version 2.1.9, serial # 0000V0000D, DEMO at 2006 04/01 07:18:33 ***

*** Example (Abaxis to LIS) sign on comment message ***

```
<STX>1H|¥&|||ABAXIS, INC. VetScan VS2^2.1.9^0000V0000D|||||D|E 1394-97|20060401071833<CR><ETX>EC<CR><LF>
<STX>2C|1|I|^Power Up Initialization Sign On|G<CR><ETX>5A<CR><LF>
<STX>3L|1|N<CR><ETX>06<CR><LF>
```

```
*** piccolo express, patient ID 22222, MALE, METLYTE 8, operator ID 1, doctor ID 1 ***
*** sample ID 44444, alternate ID 33333, birthdate 2/25/1953, MALE, phone number 789-012-3456, admission ID 88888, age 53 years
*** reference range label Patient, location 66666
```

```
<STX>1H|¥&|||ABAXIS, INC. ^piccolo xpress^2.1.9^0000V0000D| |||||D|E 1394-97|20061009072501<CR><ETX>E5<CR><LF>
<STX>2P|1|22222|44444|33333|||19530225|M|||789-012-3456|88888^53 Yrs.^|Patient| |||||66666<CR><ETX>21<CR><LF>
<STX>30|1| |^^^METLYTE 8: 4300D| |20061009072019| |1| | | |1| | | | |F<CR><ETX>40<CR><LF>
<STX>4C|1|I|^INST QC: OK      CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT:      5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>
```

```
<STX>1H|¥&|||ABAXIS, INC.^VeScan VS2^2.1.9^0000V0000D|||||D|E 1394-97|20061009072501<CR><ETX>E5<CR><LF>  
<STX>2P|122222|44444|33333|||19990630|M||||789-012-3456||88888^7 Yrs.^20060531|Dog|||||||66666<CR><ETX>21<CR><LF>  
<STX>3O|1|^|^Critical Care Profile: 0029A|[20061009072019||||2|)|||2|)|||||F<CR><ETX>42<CR><LF>  
<STX>4C|1|I|^|INST QC: OK      CHEM QC: OK|G<CR><ETX>1F<CR><LF>  
<STX>5C|2|I|^|HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>  
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>  
<STX>7C|4|I|^|Operator Renewal Required|G<CR><ETX>EA<CR><LF>
```

<STX>0C|5|I|^|^|G<CR><ETX>7D<CR><LF>
<STX>1R|1|1742-6^^LN^ALT SerP1-cCnc|20|U/L|20 to 100|N<CR><ETX>B2<CR><LF>
<STX>2R|2|2345-7^^LN^Glucose SerP1-mCnc|90|mg/dL|70 to 150|N<CR><ETX>A4<CR><LF>
<STX>3R|3|3094-0^^LN^BUN SerP1-mCnc|18|mg/dL|10 to 30|N<CR><ETX>7B<CR><LF>
<STX>4R|4|2160-0^^LN^Creat SerP1-mCnc|0.9|mg/dL|0.3 to 2.1|N<CR><ETX>0C<CR><LF>
<STX>5R|5|2951-2^^LN^Sodium SerP1-sCnc|144|mmol/L|142 to 164|N<CR><ETX>2F<CR><LF>
<STX>6R|6|2823-3^^LN^Potassium SerP1-sCnc|4.5|mmol/L|3.7 to 5.8|N<CR><ETX>83<CR><LF>
<STX>7R|7|2028-9^^LN^CO2 SerP1-sCnc|22|mmol/L|15 to 24|N<CR><ETX>ED<CR><LF>
<STX>00|2|^|^|^* QUALITY CONTROL REPORT *: 4300D||20061009072019||||1|||||1|||||||F<CR><ETX>BE<CR><LF>
<STX>1C|1|I|^|^CHEMISTRY QC: 0|G<CR><ETX>29<CR><LF>
<STX>2C|2|I|^|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>14<CR><LF>
<STX>3R|1|^|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CC<CR><LF>
<STX>4R|2|^|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CE<CR><LF>
<STX>5R|3|^|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D2<CR><LF>
<STX>6R|4|^|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>CD<CR><LF>
<STX>7R|5|^|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|6|^|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D7<CR><LF>
<STX>1R|7|^|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DD<CR><LF>
<STX>2R|8|^|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DF<CR><LF>
<STX>3L|1|N<CR><ETX>09<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID 2, matching patient ID exists so query results sent below (FEMALE sample verification issue) *****

<STX>1H|¥&||Heisenberg|Union City|510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|2<CR><ETX>EE<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401142451<CR><ETX>E7<CR><LF>
<STX>2P|1|2|44444|33333|||19530225|F|||789-012-3456||88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|^|^METLYTE 8: SS BT||20060401133051||||2|||||2|||||||Q<CR><ETX>98<CR><LF>
<STX>4C|1|I|^|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|^HEM: 181 LIP: 301 ICT: 12|G<CR><ETX>34<CR><LF>
<STX>6C|3|I|^|^@ CONFIRM LOW RECOVERIES @|G<CR><ETX>EF<CR><LF>
<STX>7C|4|I|^|^Operator Renewal Required|G<CR><ETX>6A<CR><LF>
<STX>0C|5|I|^|^|G<CR><ETX>7D<CR><LF>
<STX>1R|1|2345-7^^LN^Glucose SerP1-mCnc|20@|mg/dl|73 to 118|L<CR><ETX>06<CR><LF>
<STX>2C|1|I|419E 6666^0000 0200^@ CONFIRM LOW RECOVERIES @|I<CR><ETX>66<CR><LF>
<STX>3R|2|3094-0^^LN^BUN SerP1-mCnc|3@|mg/dl|7 to 22|L<CR><ETX>77<CR><LF>
<STX>4C|1|I|402C CCDD^0000 0200^@ CONFIRM LOW RECOVERIES @|I<CR><ETX>93<CR><LF>
<STX>5R|3|2160-0^^LN^Creat SerP1-mCnc|<0.2@|mg/dl|0.6 to 1.2|<<CR><ETX>90<CR><LF>
<STX>6C|1|I|3E29 FBE7^0000 0200^@ CONFIRM LOW RECOVERIES @|I<CR><ETX>96<CR><LF>
<STX>7R|4|2157-6^^LN^CK SerP1-cCnc|<5@|units/l|30 to 190|<<CR><ETX>AC<CR><LF>
<STX>0C|1|I|4080 0000^0000 0200^@ CONFIRM LOW RECOVERIES @|I<CR><ETX>3D<CR><LF>
<STX>1R|5|2951-2^^LN^Sodium SerP1-sCnc|^^~@|mmol/l|128 to 145|<<CR><ETX>63<CR><LF>
<STX>2C|1|I|41DE 6666^0000 0200^@ CONFIRM LOW RECOVERIES @|I<CR><ETX>71<CR><LF>
<STX>3R|6|2823-3^^LN^Potassium SerP1-sCnc|<1.5@|mmol/l|3.6 to 5.1|<<CR><ETX>FD<CR><LF>
<STX>4C|1|I|3F54 7AE1^0000 0200^@ CONFIRM LOW RECOVERIES @|I<CR><ETX>7D<CR><LF>
<STX>5R|7|2075-0^^LN^Chloride SerP1-sCnc|^^~@|mmol/l|98 to 108|<<CR><ETX>EA<CR><LF>
<STX>6C|1|I|41A8 0000^0000 0200^@ CONFIRM LOW RECOVERIES @|I<CR><ETX>4D<CR><LF>
<STX>7R|8|2028-9^^LN^CO2 SerP1-sCnc|5@|mmol/l|18 to 33|L<CR><ETX>1E<CR><LF>

*** Example (LIS to Abaxis) Query message ***

*** patient ID 3, matching patient ID exists so query results sent below ***

```
<STX>1H|¥&|||ABAXIS, INC. piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401142841<CR><ETX>EA<CR><LF>
<STX>2P|1|3|44444|33333|||19530225|M|||789-012-3456||88888^53 Yrs.^|Patient| |||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^BASIC METABOLIC PANEL: 4400D|20060401133822|||3| ||||3| |||||Q<CR><ETX>76<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^Unauthorized Operator|G<CR><ETX>8F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>7D<CR><LF>
<STX>1R|1|2345-7^LN^glucose SerPl-mCnc|96|mg/dl|73 to 118|N<CR><ETX>D4<CR><LF>
<STX>2R|2|3094-0^LN^BUN SerPl-mCnc|15|mg/dl|7 to 22|N<CR><ETX>72<CR><LF>
<STX>3R|3|17861-6^LN^Calcium SerPl-mCnc|HEM|mg/dl|8.0 to 10.3|<CR><ETX>65<CR><LF>
<STX>4C|1|I|4112 6666^OE00 0000^Hemolysis Limit Error|I<CR><ETX>6E<CR><LF>
<STX>5R|4|2160-0^LN^Creat SerPl-mCnc|0.9|mg/dl|0.6 to 1.2|N<CR><ETX>2D<CR><LF>
<STX>6R|5|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l|128 to 145|<CR><ETX>49<CR><LF>
<STX>7C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5B<CR><LF>
<STX>0R|6|2823-3^LN^Potassium SerPl-sCnc|4.0|mmol/l|3.6 to 5.1|N<CR><ETX>95<CR><LF>
<STX>1R|7|2075-0^LN^Chloride SerPl-sCnc|~|mmol/l|98 to 108|<CR><ETX>71<CR><LF>
<STX>2C|1|I|42CD 22B7^0000 0001^Chemistry Specific Error|I<CR><ETX>AC<CR><LF>
<STX>3R|8|2028-9^LN^CO2 SerPl-sCnc|26|mmol/l|18 to 33|N<CR><ETX>0E<CR><LF>
```



```

<STX>40|2|||^^^* QUALITY CONTROL REPORT *: 4400D||20060401133822|||3|||3|||3|||3|||Q<CR><ETX>CD<CR><LF>
<STX>5C|1|I|^CHEMISTRY QC:          0|G<CR><ETX>28<CR><LF>
<STX>6C|2|I|^ACCEPTABLE MINIMUM:    50|G<CR><ETX>13<CR><LF>
<STX>7R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CB<CR><LF>
<STX>0R|2|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CD<CR><LF>
<STX>1R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>2R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D4<CR><LF>
<STX>3R|5|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D0<CR><LF>
<STX>4R|6|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D6<CR><LF>
<STX>5R|7|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DC<CR><LF>
<STX>6R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DE<CR><LF>
<STX>7R|9|^LEVEL 2: PRE|100|95 to 105|<CR><ETX>94<CR><LF>
<STX>0R|10|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>48<CR><LF>
<STX>1R|11|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4D<CR><LF>
<STX>2R|12|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>57<CR><LF>
<STX>3R|13|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>45<CR><LF>
<STX>4R|14|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4D<CR><LF>
<STX>5R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4D<CR><LF>
<STX>6R|16|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|17|^LEVEL 2: 630 nm|100|95 to 105|<CR><ETX>70<CR><LF>
<STX>0L|1|N<CR><ETX>08<CR><LF>

```

*** Assume proceed with cleared database ***

*** Example (Abaxis to LIS) rotor results message ***

*** patient ID LEVEL V, CONTROL, ELECTROLYTE PANEL, operator ID 4, doctor ID 4 ***

```

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401143301<CR><ETX>E2<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^Control|||||<CR><ETX>21<CR><LF>
<STX>30|1|^ELECTROLYTE PANEL: 4200D||20060401142819||||4||||4|||||F<CR><ETX>AB<CR><LF>
<STX>4C|1|I|^INST QC: OK    CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT:    5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>2C<CR><LF>
<STX>0C|5|I|^Control Renewal Required|G<CR><ETX>7D<CR><LF>
<STX>1R|1|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F0<CR><LF>
<STX>2C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5E<CR><LF>
<STX>3R|2|2823-3^LN^Potassium SerPl-sCnc|LIP|mmol/l||<CR><ETX>43<CR><LF>
<STX>4C|1|I|4096 6666^0600 0000^Lipemic Limit Error|I<CR><ETX>70<CR><LF>
<STX>5R|3|2075-0^LN^Chloride SerPl-sCnc|^^^|mmol/l||><CR><ETX>7A<CR><LF>
<STX>6C|1|I|47C3 5000^0008 0000^Concentration > High System Limit|I<CR><ETX>7C<CR><LF>
<STX>7R|4|2028-9^LN^CO2 SerPl-sCnc|^^^|mmol/l||<CR><ETX>E0<CR><LF>
<STX>0C|1|I|41CC 0000^0020 0080^False Rate (Unexpected Reaction Occurred)|I<CR><ETX>35<CR><LF>
<STX>10|2|^* QUALITY CONTROL REPORT *: 4200D||20060401142819||||4||||4|||||F<CR><ETX>CD<CR><LF>
<STX>2C|1|I|^CHEMISTRY QC:          0|G<CR><ETX>2D<CR><LF>
<STX>3C|2|I|^ACCEPTABLE MINIMUM:    50|G<CR><ETX>10<CR><LF>
<STX>4R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>C8<CR><LF>

```

```

<STX>5R|2|^^^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CA<CR><LF>
<STX>6R|3|^^^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>CE<CR><LF>
<STX>7R|4|^^^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|5|^^^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D5<CR><LF>
<STX>1R|6|^^^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>DB<CR><LF>
<STX>2R|7|^^^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>E1<CR><LF>
<STX>3R|8|^^^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DB<CR><LF>
<STX>4R|9|^^^LEVEL 2: PRE|101|95 to 105|<CR><ETX>92<CR><LF>
<STX>5R|10|^^^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>45<CR><LF>
<STX>6R|11|^^^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|12|^^^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>54<CR><LF>
<STX>0R|13|^^^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>1R|14|^^^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>52<CR><LF>
<STX>2R|15|^^^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>52<CR><LF>
<STX>3R|16|^^^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>47<CR><LF>
<STX>4R|17|^^^LEVEL 2: 630 nm|101|95 to 105|<CR><ETX>6E<CR><LF>
<STX>5L|1|N<CR><ETX>05<CR><LF>

```

*** Example (Abaxis to LIS) rotor results message ***

*** patient ID 5, MALE, GENERAL CHEMISTRY 12, operator ID 5, doctor ID 5 ***

```

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401144000<CR><ETX>DF<CR><LF>
<STX>2P|1|5|44444|33333|||19530225|M|||789-012-3456|88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^^^GENERAL CHEMISTRY 12: 1700D|20060401143518||||5|||||5|||||||F<CR><ETX>29<CR><LF>
<STX>4C|1|I|^|^INST QC: OK    CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|^HEM: 24 LIP: 166 ICT: 1|G<CR><ETX>17<CR><LF>
<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|1751-7^^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>
<STX>2C|1|I|408C CCCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>
<STX>3R|2|6768-6^^LN^ALP SerPl-cCnc|55|units/l|53 to 128|N<CR><ETX>D1<CR><LF>
<STX>4R|3|1742-6^^LN^ALT SerPl-cCnc|HEM|units/l|10 to 47|<CR><ETX>B5<CR><LF>
<STX>5C|1|I|41E4 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>6D<CR><LF>
<STX>6R|4|1798-8^^LN^Amylase SerPl-cCnc|56|units/l|14 to 97|N<CR><ETX>98<CR><LF>
<STX>7R|5|1920-8^^LN^AST SerPl-cCnc|LIP|units/l|11 to 38|<CR><ETX>CD<CR><LF>
<STX>0C|1|I|41C4 0000^0600 0100^Lipemic Limit Error|I<CR><ETX>66<CR><LF>
<STX>1R|6|1975-2^^LN^Bilirub SerPl-mCnc|0.9|mg/dl|0.2 to 1.6|N<CR><ETX>1C<CR><LF>
<STX>2R|7|3094-0^^LN^BUN SerPl-mCnc|15|mg/dl|7 to 22|N<CR><ETX>77<CR><LF>
<STX>3R|8|17861-6^^LN^Calcium SerPl-mCnc|ICT|mg/dl|8.0 to 10.3|<CR><ETX>70<CR><LF>
<STX>4C|1|I|4112 6666^0200 0002^Icteric Limit Error|I<CR><ETX>63<CR><LF>
<STX>5R|9|2093-3^^LN^Cholest SerPl-mCnc|150|mg/dl|100 to 200|N<CR><ETX>1A<CR><LF>
<STX>6R|10|2160-0^^LN^Creat SerPl-mCnc|0.9|mg/dl|0.6 to 1.2|N<CR><ETX>5B<CR><LF>
<STX>7R|11|2345-7^^LN^Glucose SerPl-mCnc|96|mg/dl|73 to 118|N<CR><ETX>03<CR><LF>
<STX>0R|12|2885-2^^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>C3<CR><LF>
<STX>10|2|||^|^* QUALITY CONTROL REPORT *: 1700D|20060401143518||||5|||||5|||||||F<CR><ETX>CE<CR><LF>
<STX>2C|1|I|^|^CHEMISTRY QC: 0|G<CR><ETX>2D<CR><LF>

```

```

<STX>3C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>10<CR><LF>
<STX>4R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>C8<CR><LF>
<STX>5R|2|^LEVEL 1: IQC 2|102|90 to 110|<CR><ETX>CB<CR><LF>
<STX>6R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>CE<CR><LF>
<STX>7R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|5|^LEVEL 1: IQC 5|101|90 to 110|<CR><ETX>D3<CR><LF>
<STX>1R|6|^LEVEL 1: IQC 6|101|90 to 110|<CR><ETX>D6<CR><LF>
<STX>2R|7|^LEVEL 1: IQC 7|106|90 to 110|<CR><ETX>DE<CR><LF>
<STX>3R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DB<CR><LF>
<STX>4R|9|^LEVEL 2: PRE|102|95 to 105|<CR><ETX>93<CR><LF>
<STX>5R|10|^LEVEL 2: 340 nm|99|95 to 105|<CR><ETX>46<CR><LF>
<STX>6R|11|^LEVEL 2: 405 nm|100|95 to 105|<CR><ETX>69<CR><LF>
<STX>7R|12|^LEVEL 2: 467 nm|100|95 to 105|<CR><ETX>73<CR><LF>
<STX>0R|13|^LEVEL 2: 500 nm|100|95 to 105|<CR><ETX>69<CR><LF>
<STX>1R|14|^LEVEL 2: 515 nm|100|95 to 105|<CR><ETX>71<CR><LF>
<STX>2R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>52<CR><LF>
<STX>3R|16|^LEVEL 2: 600 nm|98|95 to 105|<CR><ETX>48<CR><LF>
<STX>4R|17|^LEVEL 2: 630 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>5L|1|N<CR><ETX>05<CR><LF>

```

*** Example (Abaxis to LIS) rotor results message ***

*** patient ID 6, FEMALE, LIVER PANEL PLUS, operator ID 6, doctor ID 6 ***

```

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401144524<CR><ETX>EA<CR><LF>
<STX>2P|1|6|44444|33333|||19530225|F|||789-012-3456||88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|^LIVER PANEL PLUS: 1100D|20060401144042||||6|||||6|||||||F<CR><ETX>3B<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|1751-7^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>
<STX>2C|1|I|408C CCCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>
<STX>3R|2|6768-6^LN^ALP SerPl-cCnc|55|units/l|42 to 141|N<CR><ETX>CA<CR><LF>
<STX>4R|3|1742-6^LN^ALT SerPl-cCnc|LIP|units/l|10 to 47|<CR><ETX>C0<CR><LF>
<STX>5C|1|I|41E4 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>64<CR><LF>
<STX>6R|4|1798-8^LN^Amylase SerPl-cCnc|56|units/l|14 to 97|N<CR><ETX>98<CR><LF>
<STX>7R|5|1920-8^LN^AST SerPl-cCnc|<5*|units/l|11 to 38|<CR><ETX>BF<CR><LF>
<STX>0C|1|I|41C4 0000^0001 0000^Concentration < Low Dynamic Limit|I<CR><ETX>5D<CR><LF>
<STX>1R|6|1975-2^LN^Bilirub SerPl-mCnc|0.9|mg/dl|0.2 to 1.6|N<CR><ETX>1C<CR><LF>
<STX>2R|7|2324-2^LN^GGT SerPl-cCnc|35|units/l|5 to 65|N<CR><ETX>69<CR><LF>
<STX>3R|8|2885-2^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>93<CR><LF>
<STX>40|2|^* QUALITY CONTROL REPORT *: 1100D|20060401144042||||6|||||6|||||||F<CR><ETX>BE<CR><LF>
<STX>5C|1|I|^CHEMISTRY QC: 0|G<CR><ETX>28<CR><LF>
<STX>6C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>13<CR><LF>
<STX>7R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CB<CR><LF>
<STX>0R|2|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CD<CR><LF>

```

```

<STX>1R|3|^^^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>2R|4|^^^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D4<CR><LF>
<STX>3R|5|^^^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D0<CR><LF>
<STX>4R|6|^^^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D6<CR><LF>
<STX>5R|7|^^^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DC<CR><LF>
<STX>6R|8|^^^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DE<CR><LF>
<STX>7R|9|^^^LEVEL 2: PRE|103|95 to 105|<CR><ETX>97<CR><LF>
<STX>0R|10|^^^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>48<CR><LF>
<STX>1R|11|^^^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4D<CR><LF>
<STX>2R|12|^^^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>57<CR><LF>
<STX>3R|13|^^^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>45<CR><LF>
<STX>4R|14|^^^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4D<CR><LF>
<STX>5R|15|^^^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4D<CR><LF>
<STX>6R|16|^^^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|17|^^^LEVEL 2: 630 nm|103|95 to 105|<CR><ETX>73<CR><LF>
<STX>0L|1|N<CR><ETX>08<CR><LF>

```

***** Example (Abaxis to LIS) rotor results message *****

***** patient ID LEVEL V, CONTROL, LIPID PANEL, operator ID 7, doctor ID 7 *****

```

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401145111<CR><ETX>E3<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1|^^^LIPID PANEL: 4600D||20060401144629||||7|||||7|||||||F<CR><ETX>DC<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2093-3^LN^Cholest SerPl-mCnc|ICT|mg/dl||<CR><ETX>CC<CR><LF>
<STX>2C|1|I|4353 0000^0200 0000^Icteric Limit Error|I<CR><ETX>56<CR><LF>
<STX>3R|2|2085-9^LN^HDLc SerPl-mCnc|LIP|mg/dl||<CR><ETX>3C<CR><LF>
<STX>4C|1|I|4210 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>4C<CR><LF>
<STX>5R|3|2571-8^LN^Trigl SerPl-mCnc|HEM|mg/dl||<CR><ETX>FA<CR><LF>
<STX>6C|1|I|42E6 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>71<CR><LF>
<STX>7R|4|9830-1^LN^Cholest/HDLc SerPl-mRto|^^^|||<CR><ETX>23<CR><LF>
<STX>0C|1|I|40BC CCCD^0000 0001^Chemistry Specific Error|I<CR><ETX>D6<CR><LF>
<STX>1R|5|13457-7^LN^LDLc SerPl Calc-mCnc|152c|mg/dl||<CR><ETX>25<CR><LF>
<STX>2R|6|13458-5^LN^VLDLc SerPl Calc-mCnc|23c|mg/dl||<CR><ETX>49<CR><LF>
<STX>30|2|^^^* QUALITY CONTROL REPORT *: 4600D||20060401144629||||7|||||7|||||||F<CR><ETX>D2<CR><LF>
<STX>4C|1|I|^CHEMISTRY QC: 99|G<CR><ETX>49<CR><LF>
<STX>5C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>12<CR><LF>
<STX>6R|1|^^^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CA<CR><LF>
<STX>7R|2|^^^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CC<CR><LF>
<STX>0R|3|^^^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D0<CR><LF>
<STX>1R|4|^^^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D3<CR><LF>
<STX>2R|5|^^^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D7<CR><LF>

```

<STX>3R|6|^^^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>D5<CR><LF>
<STX>4R|7|^^^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>DB<CR><LF>
<STX>5R|8|^^^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>DD<CR><LF>
<STX>6R|9|^^^LEVEL 2: PRE|103||95 to 105|<CR><ETX>96<CR><LF>
<STX>7R|10|^^^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>47<CR><LF>
<STX>0R|11|^^^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>4C<CR><LF>
<STX>1R|12|^^^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>56<CR><LF>
<STX>2R|13|^^^LEVEL 2: 500 nm|99||95 to 105|<CR><ETX>4C<CR><LF>
<STX>3R|14|^^^LEVEL 2: 515 nm|99||95 to 105|<CR><ETX>4C<CR><LF>
<STX>4R|15|^^^LEVEL 2: 550 nm|98||95 to 105|<CR><ETX>4C<CR><LF>
<STX>5R|16|^^^LEVEL 2: 600 nm|97||95 to 105|<CR><ETX>49<CR><LF>
<STX>6R|17|^^^LEVEL 2: 630 nm|103||95 to 105|<CR><ETX>72<CR><LF>
<STX>7L|1|N<CR><ETX>07<CR><LF>

*** Example (Abaxis to LIS) rotor results message ***

*** patient ID 8, HEPATIC FUNCTION PANEL, operator ID 8, doctor ID 8 ***

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401145704<CR><ETX>EB<CR><LF>
<STX>2P|1|8|44444|33333|||19530225|U|||789-012-3456|88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^^^HEPATIC FUNCTION PANEL: 4700D|20060401145222|||8|||||8|||||||F<CR><ETX>E7<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|1751-7^^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>
<STX>2C|1|I|408C CCCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>
<STX>3R|2|6768-6^^LN^ALP SerPl-cCnc|55|units/l|42 to 141|N<CR><ETX>CA<CR><LF>
<STX>4R|3|1742-6^^LN^ALT SerPl-cCnc|ICT|units/l|10 to 47|<CR><ETX>BB<CR><LF>
<STX>5C|1|I|41E8 0000^0200 0000^Icteric Limit Error|I<CR><ETX>64<CR><LF>
<STX>6R|4|1920-8^^LN^AST SerPl-cCnc|^^^|units/l|11 to 38|<CR><ETX>60<CR><LF>
<STX>7C|1|I|41C8 0000^0000 0108^Chemistry Specific Error|I<CR><ETX>87<CR><LF>
<STX>0R|5|1968-7^^LN^Bilirub Direct SerPl-mCnc|0.2|mg/dl|0.0 to 0.3|N<CR><ETX>8F<CR><LF>
<STX>1R|6|1975-2^^LN^Bilirub SerPl-mCnc|^^^|mg/dl|0.2 to 1.6|<<CR><ETX>ED<CR><LF>
<STX>2C|1|I|3F66 6666^0004 0000^Concentration < Low System Limit|I<CR><ETX>43<CR><LF>
<STX>3R|7|2885-2^^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>92<CR><LF>
<STX>40|2|||^^^* QUALITY CONTROL REPORT *: 4700D|20060401145222|||8|||||8|||||||F<CR><ETX>CC<CR><LF>
<STX>5C|1|I|^CHEMISTRY QC: 99|G<CR><ETX>4A<CR><LF>
<STX>6C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>13<CR><LF>
<STX>7R|1|^^^LEVEL 1: IQC 1|102||90 to 110|<CR><ETX>CB<CR><LF>
<STX>0R|2|^^^LEVEL 1: IQC 2|101||90 to 110|<CR><ETX>CD<CR><LF>
<STX>1R|3|^^^LEVEL 1: IQC 3|102||90 to 110|<CR><ETX>D1<CR><LF>
<STX>2R|4|^^^LEVEL 1: IQC 4|102||90 to 110|<CR><ETX>D4<CR><LF>
<STX>3R|5|^^^LEVEL 1: IQC 5|103||90 to 110|<CR><ETX>D0<CR><LF>
<STX>4R|6|^^^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>D6<CR><LF>
<STX>5R|7|^^^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>DC<CR><LF>
<STX>6R|8|^^^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>DE<CR><LF>

```

<STX>7R|9|^^^LEVEL 2: PRE|103||95 to 105|<CR><ETX>97<CR><LF>
<STX>0R|10|^^^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>48<CR><LF>
<STX>1R|11|^^^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>4D<CR><LF>
<STX>2R|12|^^^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>57<CR><LF>
<STX>3R|13|^^^LEVEL 2: 500 nm|99||95 to 105|<CR><ETX>45<CR><LF>
<STX>4R|14|^^^LEVEL 2: 515 nm|99||95 to 105|<CR><ETX>4D<CR><LF>
<STX>5R|15|^^^LEVEL 2: 550 nm|98||95 to 105|<CR><ETX>4D<CR><LF>
<STX>6R|16|^^^LEVEL 2: 600 nm|97||95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|17|^^^LEVEL 2: 630 nm|103||95 to 105|<CR><ETX>73<CR><LF>
<STX>0L|1|N<CR><ETX>08<CR><LF>

```

***** Example (Abaxis to LIS) rotor results message *****

***** patient ID 1, MALE, RENAL FUNCTION PANEL, operator ID 9, doctor ID 9 *****

```

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401150355<CR><ETX>E9<CR><LF>
<STX>2P|1|1|44444|33333|||19530225|M|||789-012-3456|88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^^^RENAL FUNCTION PANEL: 4900D||20060401145913|||9|||||9|||||||F<CR><ETX>66<CR><LF>
<STX>4C|1|I|^INST QC: OK    CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT:    5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2345-7^^LN^Glucose SerPl-mCnc|LIP|mg/dl|73 to 118|<CR><ETX>FC<CR><LF>
<STX>2C|1|I|42C0 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>64<CR><LF>
<STX>3R|2|3094-0^^LN^BUN SerPl-mCnc|15|mg/dl|7 to 22|N<CR><ETX>6B<CR><LF>
<STX>4R|3|17861-6^^LN^Calcium SerPl-mCnc|9.2|mg/dl|8.0 to 10.3|N<CR><ETX>73<CR><LF>
<STX>5R|4|2160-0^^LN^Creat SerPl-mCnc|0.9|mg/dl|0.6 to 1.2|N<CR><ETX>2D<CR><LF>
<STX>6R|5|1751-7^^LN^Alb SerPl-mCnc|4.4|g/dl|3.3 to 5.5|N<CR><ETX>F4<CR><LF>
<STX>7R|6|2777-1^^LN^Phosphate SerPl-mCnc|3.2|mg/dl|2.2 to 4.1|N<CR><ETX>F9<CR><LF>
<STX>0R|7|2951-2^^LN^Sodium SerPl-sCnc|<110*|mmol/l|128 to 145|<<CR><ETX>A1<CR><LF>
<STX>1C|1|I|4309 0000^0001 0000^Concentration < Low Dynamic Limit|I<CR><ETX>52<CR><LF>
<STX>2R|8|2823-3^^LN^Potassium SerPl-sCnc|^^^|mmol/l|3.6 to 5.1|<<CR><ETX>6F<CR><LF>
<STX>3C|1|I|4080 0000^0004 0000^Concentration < Low System Limit|I<CR><ETX>0B<CR><LF>
<STX>4R|9|2075-0^^LN^Chloride SerPl-sCnc|>135*|mmol/l|98 to 108|><CR><ETX>33<CR><LF>
<STX>5C|1|I|42CE 0000^0002 0000^Concentration > High Dynamic Limit|I<CR><ETX>BD<CR><LF>
<STX>6R|10|2028-9^^LN^CO2 SerPl-sCnc|^^^|mmol/l|18 to 33|><CR><ETX>3C<CR><LF>
<STX>7C|1|I|41D0 0000^0008 0000^Concentration > High System Limit|I<CR><ETX>70<CR><LF>
<STX>00|2|||^^^* QUALITY CONTROL REPORT *: 4900D||20060401145913|||9|||||9|||||||F<CR><ETX>DB<CR><LF>
<STX>1C|1|I|^CHEMISTRY QC:    99|G<CR><ETX>4E<CR><LF>
<STX>2C|2|I|^ACCEPTABLE MINIMUM:    50|G<CR><ETX>17<CR><LF>
<STX>3R|1|^^^LEVEL 1: IQC 1|102||90 to 110|<CR><ETX>C7<CR><LF>
<STX>4R|2|^^^LEVEL 1: IQC 2|101||90 to 110|<CR><ETX>C9<CR><LF>
<STX>5R|3|^^^LEVEL 1: IQC 3|102||90 to 110|<CR><ETX>CD<CR><LF>
<STX>6R|4|^^^LEVEL 1: IQC 4|102||90 to 110|<CR><ETX>D0<CR><LF>
<STX>7R|5|^^^LEVEL 1: IQC 5|103||90 to 110|<CR><ETX>D4<CR><LF>
<STX>0R|6|^^^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>DA<CR><LF>
<STX>1R|7|^^^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>E0<CR><LF>

```

<STX>2R|8|^^^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>E2<CR><LF>
<STX>3R|9|^^^LEVEL 2: PRE|103||95 to 105|<CR><ETX>93<CR><LF>
<STX>4R|10|^^^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>44<CR><LF>
<STX>5R|11|^^^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>49<CR><LF>
<STX>6R|12|^^^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>53<CR><LF>
<STX>7R|13|^^^LEVEL 2: 500 nm|99||95 to 105|<CR><ETX>49<CR><LF>
<STX>0R|14|^^^LEVEL 2: 515 nm|99||95 to 105|<CR><ETX>51<CR><LF>
<STX>1R|15|^^^LEVEL 2: 550 nm|98||95 to 105|<CR><ETX>51<CR><LF>
<STX>2R|16|^^^LEVEL 2: 600 nm|97||95 to 105|<CR><ETX>4E<CR><LF>
<STX>3R|17|^^^LEVEL 2: 630 nm|103||95 to 105|<CR><ETX>6F<CR><LF>
<STX>4L|1|N<CR><ETX>04<CR><LF>

***** Example (Abaxis to LIS) rotor results message *****

***** patient ID LEVEL IX, CONTROL, METLYTE 8, operator ID 0, doctor ID 0 *****

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D||||||D|E 1394-97|20060401150945<CR><ETX>EE<CR><LF>
<STX>2P|1|LEVEL IX|||||U|||||^Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1|||^METLYTE 8: 4300D|20060401150503||||0|||||0|||||||F<CR><ETX>39<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2345-7^^LN^Glucose SerPl-mCnc|HEM|mg/dl||<CR><ETX>CA<CR><LF>
<STX>2C|1|I|42BF 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>82<CR><LF>
<STX>3R|2|3094-0^^LN^BUN SerPl-mCnc|15|mg/dl||<CR><ETX>5F<CR><LF>
<STX>4R|3|2160-0^^LN^Creat SerPl-mCnc|LIP|mg/dl||<CR><ETX>E3<CR><LF>
<STX>5C|1|I|3F66 6666^0600 0000^Lipemic Limit Error|I<CR><ETX>83<CR><LF>
<STX>6R|4|2157-6^^LN^CK SerPl-cCnc|115|units/l||<CR><ETX>34<CR><LF>
<STX>7R|5|2951-2^^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F2<CR><LF>
<STX>0C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5C<CR><LF>
<STX>1R|6|2823-3^^LN^Potassium SerPl-sCnc|4.0|mmol/l||<CR><ETX>FA<CR><LF>
<STX>2R|7|2075-0^^LN^Chloride SerPl-sCnc|^|mmol/l||><CR><ETX>83<CR><LF>
<STX>3C|1|I|461C 4000^0008 0000^Concentration > High System Limit|I<CR><ETX>75<CR><LF>
<STX>4R|8|2028-9^^LN^C02 SerPl-sCnc|26|mmol/l||<CR><ETX>CF<CR><LF>
<STX>50|2|||^* QUALITY CONTROL REPORT *: 4300D|20060401150503||||0|||||0|||||||F<CR><ETX>B7<CR><LF>
<STX>6C|1|I|^CHEMISTRY QC: 99|G<CR><ETX>4B<CR><LF>
<STX>7C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>14<CR><LF>
<STX>0R|1|^^^LEVEL 1: IQC 1|102||90 to 110|<CR><ETX>CC<CR><LF>
<STX>1R|2|^^^LEVEL 1: IQC 2|101||90 to 110|<CR><ETX>CE<CR><LF>
<STX>2R|3|^^^LEVEL 1: IQC 3|102||90 to 110|<CR><ETX>D2<CR><LF>
<STX>3R|4|^^^LEVEL 1: IQC 4|102||90 to 110|<CR><ETX>CD<CR><LF>
<STX>4R|5|^^^LEVEL 1: IQC 5|103||90 to 110|<CR><ETX>D1<CR><LF>
<STX>5R|6|^^^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>D7<CR><LF>
<STX>6R|7|^^^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>DD<CR><LF>
<STX>7R|8|^^^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>DF<CR><LF>
<STX>0R|9|^^^LEVEL 2: PRE|100||95 to 105|<CR><ETX>95<CR><LF>

<STX>1R|10|^^^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>49<CR><LF>
<STX>2R|11|^^^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>4E<CR><LF>
<STX>3R|12|^^^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>50<CR><LF>
<STX>4R|13|^^^LEVEL 2: 500 nm|99||95 to 105|<CR><ETX>46<CR><LF>
<STX>5R|14|^^^LEVEL 2: 515 nm|99||95 to 105|<CR><ETX>4E<CR><LF>
<STX>6R|15|^^^LEVEL 2: 550 nm|98||95 to 105|<CR><ETX>4E<CR><LF>
<STX>7R|16|^^^LEVEL 2: 600 nm|97||95 to 105|<CR><ETX>4B<CR><LF>
<STX>0R|17|^^^LEVEL 2: 630 nm|100||95 to 105|<CR><ETX>71<CR><LF>
<STX>1L|1|N<CR><ETX>09<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID all (since restart), all rotors ran have matching query results dumped below *****

<STX>1H|¥^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|all<CR><ETX>F5<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D||||||D|E 1394-97|20060401143301<CR><ETX>E2<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^|^Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1|^^^ELECTROLYTE PANEL: 4200D||20060401142819||||4|||||4|||||||Q<CR><ETX>AB<CR><LF>
<STX>4C|1|I|^|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|^|G<CR><ETX>2C<CR><LF>
<STX>0C|5|I|^|^Control Renewal Required|G<CR><ETX>7D<CR><LF>
<STX>1R|1|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F0<CR><LF>
<STX>2C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5E<CR><LF>
<STX>3R|2|2823-3^LN^Potassium SerPl-sCnc|LIP|mmol/l||<CR><ETX>43<CR><LF>
<STX>4C|1|I|4096 6666^0600 0000^Lipemic Limit Error|I<CR><ETX>70<CR><LF>
<STX>5R|3|2075-0^LN^Chloride SerPl-sCnc|^|^mmol/l||><CR><ETX>7A<CR><LF>
<STX>6C|1|I|47C3 5000^0008 0000^Concentration > High System Limit|I<CR><ETX>7C<CR><LF>
<STX>7R|4|2028-9^LN^CO2 SerPl-sCnc|^|^mmol/l||<CR><ETX>E0<CR><LF>
<STX>0C|1|I|41CC 0000^0020 0080^False Rate (Unexpected Reaction Occurred)|I<CR><ETX>35<CR><LF>
<STX>10|2|^^^* QUALITY CONTROL REPORT *: 4200D||20060401142819||||4|||||4|||||||Q<CR><ETX>CD<CR><LF>
<STX>2C|1|I|^|^CHEMISTRY QC: 0|G<CR><ETX>2D<CR><LF>
<STX>3C|2|I|^|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>10<CR><LF>
<STX>4R|1|^|^LEVEL 1: IQC 1|102||90 to 110|<CR><ETX>C8<CR><LF>
<STX>5R|2|^|^LEVEL 1: IQC 2|101||90 to 110|<CR><ETX>CA<CR><LF>
<STX>6R|3|^|^LEVEL 1: IQC 3|102||90 to 110|<CR><ETX>CE<CR><LF>
<STX>7R|4|^|^LEVEL 1: IQC 4|102||90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|5|^|^LEVEL 1: IQC 5|103||90 to 110|<CR><ETX>D5<CR><LF>
<STX>1R|6|^|^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>DB<CR><LF>
<STX>2R|7|^|^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>E1<CR><LF>
<STX>3R|8|^|^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>DB<CR><LF>
<STX>4R|9|^|^LEVEL 2: PRE|101||95 to 105|<CR><ETX>92<CR><LF>
<STX>5R|10|^|^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>45<CR><LF>
<STX>6R|11|^|^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|12|^|^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>54<CR><LF>

<STX>0R|13|^^^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>1R|14|^^^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>52<CR><LF>
<STX>2R|15|^^^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>52<CR><LF>
<STX>3R|16|^^^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>47<CR><LF>
<STX>4R|17|^^^LEVEL 2: 630 nm|101|95 to 105|<CR><ETX>6E<CR><LF>
<STX>5L|1|N<CR><ETX>05<CR><LF>

<STX>1H|¥&|||ABAXIS, INC. ^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401144000<CR><ETX>DF<CR><LF>
<STX>2P|1|5|44444|33333|||19530225|M|||789-012-3456|88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^^^GENERAL CHEMISTRY 12: 1700D||20060401143518||||5||||5|||||||Q<CR><ETX>29<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: 24 LIP: 166 ICT: 1|G<CR><ETX>17<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|1751-7^^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>
<STX>2C|1|I|408C CCCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>
<STX>3R|2|6768-6^^LN^ALP SerPl-cCnc|55|units/l|53 to 128|N<CR><ETX>D1<CR><LF>
<STX>4R|3|1742-6^^LN^ALT SerPl-cCnc|HEM|units/l|10 to 47|<CR><ETX>B5<CR><LF>
<STX>5C|1|I|41E4 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>6D<CR><LF>
<STX>6R|4|1798-8^^LN^Amylase SerPl-cCnc|56|units/l|14 to 97|N<CR><ETX>98<CR><LF>
<STX>7R|5|1920-8^^LN^AST SerPl-cCnc|LIP|units/l|11 to 38|<CR><ETX>CD<CR><LF>
<STX>0C|1|I|41C4 0000^0600 0100^Lipemic Limit Error|I<CR><ETX>66<CR><LF>
<STX>1R|6|1975-2^^LN^Bilirub SerPl-mCnc|0.9|mg/dl|0.2 to 1.6|N<CR><ETX>1C<CR><LF>
<STX>2R|7|3094-0^^LN^BUN SerPl-mCnc|15|mg/dl|7 to 22|N<CR><ETX>77<CR><LF>
<STX>3R|8|17861-6^^LN^Calcium SerPl-mCnc|ICT|mg/dl|8.0 to 10.3|<CR><ETX>70<CR><LF>
<STX>4C|1|I|4112 6666^0200 0002^Icteric Limit Error|I<CR><ETX>63<CR><LF>
<STX>5R|9|2093-3^^LN^Cholest SerPl-mCnc|150|mg/dl|100 to 200|N<CR><ETX>1A<CR><LF>
<STX>6R|10|2160-0^^LN^Creat SerPl-mCnc|0.9|mg/dl|0.6 to 1.2|N<CR><ETX>5B<CR><LF>
<STX>7R|11|2345-7^^LN^Glucose SerPl-mCnc|96|mg/dl|73 to 118|N<CR><ETX>03<CR><LF>
<STX>0R|12|2885-2^^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>C3<CR><LF>
<STX>10|2|||^^^* QUALITY CONTROL REPORT *: 1700D||20060401143518||||5||||5|||||||Q<CR><ETX>CE<CR><LF>
<STX>2C|1|I|^CHEMISTRY QC: 0|G<CR><ETX>2D<CR><LF>
<STX>3C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>10<CR><LF>
<STX>4R|1|^^^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>C8<CR><LF>
<STX>5R|2|^^^LEVEL 1: IQC 2|102|90 to 110|<CR><ETX>CB<CR><LF>
<STX>6R|3|^^^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>CE<CR><LF>
<STX>7R|4|^^^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|5|^^^LEVEL 1: IQC 5|101|90 to 110|<CR><ETX>D3<CR><LF>
<STX>1R|6|^^^LEVEL 1: IQC 6|101|90 to 110|<CR><ETX>D6<CR><LF>
<STX>2R|7|^^^LEVEL 1: IQC 7|106|90 to 110|<CR><ETX>DE<CR><LF>
<STX>3R|8|^^^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DB<CR><LF>
<STX>4R|9|^^^LEVEL 2: PRE|102|95 to 105|<CR><ETX>93<CR><LF>
<STX>5R|10|^^^LEVEL 2: 340 nm|99|95 to 105|<CR><ETX>46<CR><LF>
<STX>6R|11|^^^LEVEL 2: 405 nm|100|95 to 105|<CR><ETX>69<CR><LF>
<STX>7R|12|^^^LEVEL 2: 467 nm|100|95 to 105|<CR><ETX>73<CR><LF>
<STX>0R|13|^^^LEVEL 2: 500 nm|100|95 to 105|<CR><ETX>69<CR><LF>
<STX>1R|14|^^^LEVEL 2: 515 nm|100|95 to 105|<CR><ETX>71<CR><LF>

<STX>2R|15|^^^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>52<CR><LF>
<STX>3R|16|^^^LEVEL 2: 600 nm|98|95 to 105|<CR><ETX>48<CR><LF>
<STX>4R|17|^^^LEVEL 2: 630 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>5L|1|N<CR><ETX>05<CR><LF>

<STX>1H|¥&|||ABAXIS, INC. ^piccolo xpress^2.1.9^0000V0000D||||||D|E 1394-97|20060401144524<CR><ETX>EA<CR><LF>
<STX>2P|1|6|44444|33333||19530225|F|||789-012-3456|88888^53 Yrs.^Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^LIVER PANEL PLUS: 1100D|20060401144042|||6|||||6||||||Q<CR><ETX>3B<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|1751-7^^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>
<STX>2C|1|I|408C CCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>
<STX>3R|2|6768-6^^LN^ALP SerPl-cCnc|55|units/l|42 to 141|N<CR><ETX>CA<CR><LF>
<STX>4R|3|1742-6^^LN^ALT SerPl-cCnc|LIP|units/l|10 to 47|<CR><ETX>C0<CR><LF>
<STX>5C|1|I|41E4 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>64<CR><LF>
<STX>6R|4|1798-8^^LN^Amylase SerPl-cCnc|56|units/l|14 to 97|N<CR><ETX>98<CR><LF>
<STX>7R|5|1920-8^^LN^AST SerPl-cCnc|<5*|units/l|11 to 38|<<CR><ETX>BF<CR><LF>
<STX>0C|1|I|41C4 0000^0001 0000^Concentration < Low Dynamic Limit|I<CR><ETX>5D<CR><LF>
<STX>1R|6|1975-2^^LN^Bilirub SerPl-mCnc|0.9|mg/dl|0.2 to 1.6|N<CR><ETX>1C<CR><LF>
<STX>2R|7|2324-2^^LN^GGT SerPl-cCnc|35|units/l|5 to 65|N<CR><ETX>69<CR><LF>
<STX>3R|8|2885-2^^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>93<CR><LF>
<STX>40|2|||^* QUALITY CONTROL REPORT *: 1100D|20060401144042|||6|||||6||||||Q<CR><ETX>BE<CR><LF>
<STX>5C|1|I|^CHEMISTRY QC: 0|G<CR><ETX>28<CR><LF>
<STX>6C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>13<CR><LF>
<STX>7R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CB<CR><LF>
<STX>0R|2|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CD<CR><LF>
<STX>1R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>2R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D4<CR><LF>
<STX>3R|5|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D0<CR><LF>
<STX>4R|6|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D6<CR><LF>
<STX>5R|7|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DC<CR><LF>
<STX>6R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DE<CR><LF>
<STX>7R|9|^LEVEL 2: PRE|103|95 to 105|<CR><ETX>97<CR><LF>
<STX>0R|10|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>48<CR><LF>
<STX>1R|11|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4D<CR><LF>
<STX>2R|12|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>57<CR><LF>
<STX>3R|13|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>45<CR><LF>
<STX>4R|14|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4D<CR><LF>
<STX>5R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4D<CR><LF>
<STX>6R|16|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|17|^LEVEL 2: 630 nm|103|95 to 105|<CR><ETX>73<CR><LF>
<STX>0L|1|N<CR><ETX>08<CR><LF>

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D||||||D|E 1394-97|20060401145111<CR><ETX>E3<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^|Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1|^^^LIPID PANEL: 4600D||20060401144629||||7|||||7|||||||Q<CR><ETX>DC<CR><LF>
<STX>4C|1|I|^|INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2093-3^LN^Cholest SerPl-mCnc|ICT|mg/dl||<CR><ETX>CC<CR><LF>
<STX>2C|1|I|4353 0000^0200 0000^Icteric Limit Error|I<CR><ETX>56<CR><LF>
<STX>3R|2|2085-9^LN^HDLc SerPl-mCnc|LIP|mg/dl||<CR><ETX>3C<CR><LF>
<STX>4C|1|I|4210 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>4C<CR><LF>
<STX>5R|3|2571-8^LN^Trigl SerPl-mCnc|HEM|mg/dl||<CR><ETX>FA<CR><LF>
<STX>6C|1|I|42E6 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>71<CR><LF>
<STX>7R|4|9830-1^LN^Cholest/HDLc SerPl-mRto|^^^||<CR><ETX>23<CR><LF>
<STX>0C|1|I|40BC CCCD^0000 0001^Chemistry Specific Error|I<CR><ETX>D6<CR><LF>
<STX>1R|5|13457-7^LN^LDLc SerPl Calc-mCnc|152c|mg/dl||<CR><ETX>25<CR><LF>
<STX>2R|6|13458-5^LN^VLDLc SerPl Calc-mCnc|23c|mg/dl||<CR><ETX>49<CR><LF>
<STX>30|2|^^^* QUALITY CONTROL REPORT *: 4600D||20060401144629||||7|||||7|||||||Q<CR><ETX>D2<CR><LF>
<STX>4C|1|I|^|CHEMISTRY QC: 99|G<CR><ETX>49<CR><LF>
<STX>5C|2|I|^|ACCEPTABLE MINIMUM: 50|G<CR><ETX>12<CR><LF>
<STX>6R|1|^^^LEVEL 1: IQC 1|102||90 to 110|<CR><ETX>CA<CR><LF>
<STX>7R|2|^^^LEVEL 1: IQC 2|101||90 to 110|<CR><ETX>CC<CR><LF>
<STX>0R|3|^^^LEVEL 1: IQC 3|102||90 to 110|<CR><ETX>D0<CR><LF>
<STX>1R|4|^^^LEVEL 1: IQC 4|102||90 to 110|<CR><ETX>D3<CR><LF>
<STX>2R|5|^^^LEVEL 1: IQC 5|103||90 to 110|<CR><ETX>D7<CR><LF>
<STX>3R|6|^^^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>D5<CR><LF>
<STX>4R|7|^^^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>DB<CR><LF>
<STX>5R|8|^^^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>DD<CR><LF>
<STX>6R|9|^^^LEVEL 2: PRE|103||95 to 105|<CR><ETX>96<CR><LF>
<STX>7R|10|^^^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>47<CR><LF>
<STX>0R|11|^^^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>4C<CR><LF>
<STX>1R|12|^^^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>56<CR><LF>
<STX>2R|13|^^^LEVEL 2: 500 nm|99||95 to 105|<CR><ETX>4C<CR><LF>
<STX>3R|14|^^^LEVEL 2: 515 nm|99||95 to 105|<CR><ETX>4C<CR><LF>
<STX>4R|15|^^^LEVEL 2: 550 nm|98||95 to 105|<CR><ETX>4C<CR><LF>
<STX>5R|16|^^^LEVEL 2: 600 nm|97||95 to 105|<CR><ETX>49<CR><LF>
<STX>6R|17|^^^LEVEL 2: 630 nm|103||95 to 105|<CR><ETX>72<CR><LF>
<STX>7L|1|N<CR><ETX>07<CR><LF>

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D||||||D|E 1394-97|20060401145704<CR><ETX>EB<CR><LF>
<STX>2P|1|8|44444|33333|||19530225|U|||789-012-3456|88888^53 Yrs.^Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^^^HEPATIC FUNCTION PANEL: 4700D|20060401145222||||8|||||8|||||||Q<CR><ETX>E7<CR><LF>
<STX>4C|1|I|^|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|1751-7^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>
<STX>2C|1|I|408C CCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>
<STX>3R|2|6768-6^LN^ALP SerPl-cCnc|55|units/l|42 to 141|N<CR><ETX>CA<CR><LF>
<STX>4R|3|1742-6^LN^ALT SerPl-cCnc|ICT|units/l|10 to 47|<CR><ETX>BB<CR><LF>
<STX>5C|1|I|41E8 0000^0200 0000^Icteric Limit Error|I<CR><ETX>64<CR><LF>
<STX>6R|4|1920-8^LN^AST SerPl-cCnc|^|^|units/l|11 to 38|<CR><ETX>60<CR><LF>
<STX>7C|1|I|41C8 0000^0000 0108^Chemistry Specific Error|I<CR><ETX>87<CR><LF>
<STX>0R|5|1968-7^LN^Bilirub Direct SerPl-mCnc|0.2|mg/dl|0.0 to 0.3|N<CR><ETX>8F<CR><LF>
<STX>1R|6|1975-2^LN^Bilirub SerPl-mCnc|^|^|mg/dl|0.2 to 1.6|<<CR><ETX>ED<CR><LF>
<STX>2C|1|I|3F66 6666^0004 0000^Concentration < Low System Limit|I<CR><ETX>43<CR><LF>
<STX>3R|7|2885-2^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>92<CR><LF>
<STX>40|2|||^|^* QUALITY CONTROL REPORT *: 4700D|20060401145222||||8|||||8|||||||Q<CR><ETX>CC<CR><LF>
<STX>5C|1|I|^|^CHEMISTRY QC: 99|G<CR><ETX>4A<CR><LF>
<STX>6C|2|I|^|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>13<CR><LF>
<STX>7R|1|^|^LEVEL 1: IQC 1|102||90 to 110|<CR><ETX>CB<CR><LF>
<STX>0R|2|^|^LEVEL 1: IQC 2|101||90 to 110|<CR><ETX>CD<CR><LF>
<STX>1R|3|^|^LEVEL 1: IQC 3|102||90 to 110|<CR><ETX>D1<CR><LF>
<STX>2R|4|^|^LEVEL 1: IQC 4|102||90 to 110|<CR><ETX>D4<CR><LF>
<STX>3R|5|^|^LEVEL 1: IQC 5|103||90 to 110|<CR><ETX>D0<CR><LF>
<STX>4R|6|^|^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>D6<CR><LF>
<STX>5R|7|^|^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>DC<CR><LF>
<STX>6R|8|^|^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>DE<CR><LF>
<STX>7R|9|^|^LEVEL 2: PRE|103||95 to 105|<CR><ETX>97<CR><LF>
<STX>0R|10|^|^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>48<CR><LF>
<STX>1R|11|^|^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>4D<CR><LF>
<STX>2R|12|^|^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>57<CR><LF>
<STX>3R|13|^|^LEVEL 2: 500 nm|99||95 to 105|<CR><ETX>45<CR><LF>
<STX>4R|14|^|^LEVEL 2: 515 nm|99||95 to 105|<CR><ETX>4D<CR><LF>
<STX>5R|15|^|^LEVEL 2: 550 nm|98||95 to 105|<CR><ETX>4D<CR><LF>
<STX>6R|16|^|^LEVEL 2: 600 nm|97||95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|17|^|^LEVEL 2: 630 nm|103||95 to 105|<CR><ETX>73<CR><LF>
<STX>0L|1|N<CR><ETX>08<CR><LF>

<STX>1H|¥&|||ABAXIS, INC. ^piccolo xpress^2.1.9^0000V0000D||||||D|E 1394-97|20060401150355<CR><ETX>E9<CR><LF>
<STX>2P|1|1|44444|33333|||19530225|M|||789-012-3456||88888^53 Yrs.^Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^^^RENAL FUNCTION PANEL: 4900D||20060401145913|||9|||||9|||||||Q<CR><ETX>66<CR><LF>
<STX>4C|1|I|^|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2345-7^LN^Glucose SerPl-mCnc|LIP|mg/dl|73 to 118|<CR><ETX>FC<CR><LF>
<STX>2C|1|I|42C0 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>64<CR><LF>
<STX>3R|2|3094-0^LN^BUN SerPl-mCnc|15|mg/dl|7 to 22|N<CR><ETX>6B<CR><LF>
<STX>4R|3|17861-6^LN^Calcium SerPl-mCnc|9.2|mg/dl|8.0 to 10.3|N<CR><ETX>73<CR><LF>
<STX>5R|4|2160-0^LN^Creat SerPl-mCnc|0.9|mg/dl|0.6 to 1.2|N<CR><ETX>2D<CR><LF>
<STX>6R|5|1751-7^LN^Alb SerPl-mCnc|4.4|g/dl|3.3 to 5.5|N<CR><ETX>F4<CR><LF>
<STX>7R|6|2777-1^LN^Phosphate SerPl-mCnc|3.2|mg/dl|2.2 to 4.1|N<CR><ETX>F9<CR><LF>
<STX>0R|7|2951-2^LN^Sodium SerPl-sCnc|<110*|mmol/l|128 to 145|<<CR><ETX>A1<CR><LF>
<STX>1C|1|I|4309 0000^0001 0000^Concentration < Low Dynamic Limit|I<CR><ETX>52<CR><LF>
<STX>2R|8|2823-3^LN^Potassium SerPl-sCnc|^^^|mmol/l|3.6 to 5.1|<<CR><ETX>6F<CR><LF>
<STX>3C|1|I|4080 0000^0004 0000^Concentration < Low System Limit|I<CR><ETX>0B<CR><LF>
<STX>4R|9|2075-0^LN^Chloride SerPl-sCnc|>135*|mmol/l|98 to 108|><CR><ETX>33<CR><LF>
<STX>5C|1|I|42CE 0000^0002 0000^Concentration > High Dynamic Limit|I<CR><ETX>BD<CR><LF>
<STX>6R|10|2028-9^LN^CO2 SerPl-sCnc|^^^|mmol/l|18 to 33|><CR><ETX>3C<CR><LF>
<STX>7C|1|I|41D0 0000^0008 0000^Concentration > High System Limit|I<CR><ETX>70<CR><LF>
<STX>00|2|||^|^* QUALITY CONTROL REPORT *: 4900D||20060401145913|||9|||||9|||||||Q<CR><ETX>DB<CR><LF>
<STX>1C|1|I|^|^CHEMISTRY QC: 99|G<CR><ETX>4E<CR><LF>
<STX>2C|2|I|^|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>17<CR><LF>
<STX>3R|1|^|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>C7<CR><LF>
<STX>4R|2|^|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>C9<CR><LF>
<STX>5R|3|^|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>CD<CR><LF>
<STX>6R|4|^|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D0<CR><LF>
<STX>7R|5|^|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D4<CR><LF>
<STX>0R|6|^|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>DA<CR><LF>
<STX>1R|7|^|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>E0<CR><LF>
<STX>2R|8|^|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>E2<CR><LF>
<STX>3R|9|^|^LEVEL 2: PRE|103|95 to 105|<CR><ETX>93<CR><LF>
<STX>4R|10|^|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>44<CR><LF>
<STX>5R|11|^|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>49<CR><LF>
<STX>6R|12|^|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>53<CR><LF>
<STX>7R|13|^|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>49<CR><LF>
<STX>0R|14|^|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>51<CR><LF>
<STX>1R|15|^|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>51<CR><LF>
<STX>2R|16|^|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>4E<CR><LF>
<STX>3R|17|^|^LEVEL 2: 630 nm|103|95 to 105|<CR><ETX>6F<CR><LF>
<STX>4L|1|N<CR><ETX>04<CR><LF>

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401150945<CR><ETX>EE<CR><LF>
<STX>2P|1|LEVEL 1X|||||U|||||^|Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1|||^|^METLYTE 8: 4300D|20060401150503|||0|||||0|||||||Q<CR><ETX>39<CR><LF>
<STX>4C|1|I|^|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2345-7^LN^Glucose SerPl-mCnc|HEM|mg/dl||<CR><ETX>CA<CR><LF>
<STX>2C|1|I|42BF 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>82<CR><LF>
<STX>3R|2|3094-0^LN^BUN SerPl-mCnc|15|mg/dl||<CR><ETX>5F<CR><LF>
<STX>4R|3|2160-0^LN^Creat SerPl-mCnc|LIP|mg/dl||<CR><ETX>E3<CR><LF>
<STX>5C|1|I|3F66 6666^0600 0000^Lipemic Limit Error|I<CR><ETX>83<CR><LF>
<STX>6R|4|2157-6^LN^CK SerPl-cCnc|115|units/l||<CR><ETX>34<CR><LF>
<STX>7R|5|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F2<CR><LF>
<STX>0C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5C<CR><LF>
<STX>1R|6|2823-3^LN^Potassium SerPl-sCnc|4.0|mmol/l||<CR><ETX>FA<CR><LF>
<STX>2R|7|2075-0^LN^Chloride SerPl-sCnc|^|mmol/l||><CR><ETX>83<CR><LF>
<STX>3C|1|I|461C 4000^0008 0000^Concentration > High System Limit|I<CR><ETX>75<CR><LF>
<STX>4R|8|2028-9^LN^CO2 SerPl-sCnc|26|mmol/l||<CR><ETX>CF<CR><LF>
<STX>50|2|||^|^* QUALITY CONTROL REPORT *: 4300D|20060401150503|||0|||||0|||||||Q<CR><ETX>B7<CR><LF>
<STX>6C|1|I|^|^CHEMISTRY QC: 99|G<CR><ETX>4B<CR><LF>
<STX>7C|2|I|^|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>14<CR><LF>
<STX>0R|1|^|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CC<CR><LF>
<STX>1R|2|^|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CE<CR><LF>
<STX>2R|3|^|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D2<CR><LF>
<STX>3R|4|^|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>CD<CR><LF>
<STX>4R|5|^|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D1<CR><LF>
<STX>5R|6|^|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D7<CR><LF>
<STX>6R|7|^|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DD<CR><LF>
<STX>7R|8|^|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DF<CR><LF>
<STX>0R|9|^|^LEVEL 2: PRE|100|95 to 105|<CR><ETX>95<CR><LF>
<STX>1R|10|^|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>49<CR><LF>
<STX>2R|11|^|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>3R|12|^|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>50<CR><LF>
<STX>4R|13|^|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>46<CR><LF>
<STX>5R|14|^|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>6R|15|^|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4E<CR><LF>
<STX>7R|16|^|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>4B<CR><LF>
<STX>0R|17|^|^LEVEL 2: 630 nm|100|95 to 105|<CR><ETX>71<CR><LF>
<STX>1L|1|N<CR><ETX>09<CR><LF>

*** Example (LIS to Abaxis) Query message ***

*** patient ID 1, matching query results for ID 1 dumped below ***

```
<STX>1H|¥^&|||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|1<CR><ETX>ED<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥^&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401150355<CR><ETX>E9<CR><LF>
<STX>2P|1|1|44444|33333|||19530225|M|||789-012-3456||88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|1|^RENAL FUNCTION PANEL: 4900D||20060401145913|||9|||9|||Q<CR><ETX>66<CR><LF>
<STX>4C|1|I|^INST QC: OK   CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2345-7^LN^Glucose SerPl-mCnc|LIP|mg/dl|73 to 118|<CR><ETX>FC<CR><LF>
<STX>2C|1|I|42C0 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>64<CR><LF>
<STX>3R|2|3094-0^LN^BUN SerPl-mCnc|15|mg/dl|7 to 22|N<CR><ETX>6B<CR><LF>
<STX>4R|3|17861-6^LN^Calcium SerPl-mCnc|9.2|mg/dl|8.0 to 10.3|N<CR><ETX>73<CR><LF>
<STX>5R|4|2160-0^LN^Creat SerPl-mCnc|0.9|mg/dl|0.6 to 1.2|N<CR><ETX>2D<CR><LF>
<STX>6R|5|1751-7^LN^Alb SerPl-mCnc|4.4|g/dl|3.3 to 5.5|N<CR><ETX>F4<CR><LF>
<STX>7R|6|2777-1^LN^Phosphate SerPl-mCnc|3.2|mg/dl|2.2 to 4.1|N<CR><ETX>F9<CR><LF>
<STX>0R|7|2951-2^LN^Sodium SerPl-sCnc|<110*|mmol/l|128 to 145|<<CR><ETX>A1<CR><LF>
<STX>1C|1|I|4309 0000^0001 0000^Concentration < Low Dynamic Limit|I<CR><ETX>52<CR><LF>
<STX>2R|8|2823-3^LN^Potassium SerPl-sCnc|^^^|mmol/l|3.6 to 5.1|<<CR><ETX>6F<CR><LF>
<STX>3C|1|I|4080 0000^0004 0000^Concentration < Low System Limit|I<CR><ETX>0B<CR><LF>
<STX>4R|9|2075-0^LN^Chloride SerPl-sCnc|>135*|mmol/l|98 to 108|><CR><ETX>33<CR><LF>
<STX>5C|1|I|42CE 0000^0002 0000^Concentration > High Dynamic Limit|I<CR><ETX>BD<CR><LF>
<STX>6R|10|2028-9^LN^CO2 SerPl-sCnc|^^^|mmol/l|18 to 33|><CR><ETX>3C<CR><LF>
<STX>7C|1|I|41D0 0000^0008 0000^Concentration > High System Limit|I<CR><ETX>70<CR><LF>
<STX>00|2|1|^* QUALITY CONTROL REPORT *: 4900D||20060401145913|||9|||9|||Q<CR><ETX>DB<CR><LF>
<STX>1C|1|I|^CHEMISTRY QC: 99|G<CR><ETX>4E<CR><LF>
<STX>2C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>17<CR><LF>
<STX>3R|1|^LEVEL 1: IQC 1|102||90 to 110|<CR><ETX>C7<CR><LF>
<STX>4R|2|^LEVEL 1: IQC 2|101||90 to 110|<CR><ETX>C9<CR><LF>
<STX>5R|3|^LEVEL 1: IQC 3|102||90 to 110|<CR><ETX>CD<CR><LF>
<STX>6R|4|^LEVEL 1: IQC 4|102||90 to 110|<CR><ETX>D0<CR><LF>
<STX>7R|5|^LEVEL 1: IQC 5|103||90 to 110|<CR><ETX>D4<CR><LF>
<STX>0R|6|^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>DA<CR><LF>
<STX>1R|7|^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>E0<CR><LF>
<STX>2R|8|^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>E2<CR><LF>
<STX>3R|9|^LEVEL 2: PRE|103||95 to 105|<CR><ETX>93<CR><LF>
<STX>4R|10|^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>44<CR><LF>
<STX>5R|11|^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>49<CR><LF>
<STX>6R|12|^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>53<CR><LF>
<STX>7R|13|^LEVEL 2: 500 nm|99||95 to 105|<CR><ETX>49<CR><LF>
<STX>0R|14|^LEVEL 2: 515 nm|99||95 to 105|<CR><ETX>51<CR><LF>
<STX>1R|15|^LEVEL 2: 550 nm|98||95 to 105|<CR><ETX>51<CR><LF>
<STX>2R|16|^LEVEL 2: 600 nm|97||95 to 105|<CR><ETX>4E<CR><LF>
<STX>3R|17|^LEVEL 2: 630 nm|103||95 to 105|<CR><ETX>6F<CR><LF>
<STX>4L|1|N<CR><ETX>04<CR><LF>
```

*** Example (LIS to Abaxis) Query message ***

*** patient ID level v, matching query results for LEVEL V CONTROL rotors dumped below ***

```
<STX>1H|¥&|||Heisenberg|Union City|510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|level v<CR><ETX>6A<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401143301<CR><ETX>E2<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^|Control|||||||<CR><ETX>21<CR><LF>
<STX>3O|1|||^ELECTROLYTE PANEL: 4200D|20060401142819||||4||||4|||||||Q<CR><ETX>AB<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>2C<CR><LF>
<STX>0C|5|I|^Control Renewal Required|G<CR><ETX>7D<CR><LF>
<STX>1R|1|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F0<CR><LF>
<STX>2C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5E<CR><LF>
<STX>3R|2|2823-3^LN^Potassium SerPl-sCnc|LIP|mmol/l||<CR><ETX>43<CR><LF>
<STX>4C|1|I|4096 6666^0600 0000^Lipemic Limit Error|I<CR><ETX>70<CR><LF>
<STX>5R|3|2075-0^LN^Chloride SerPl-sCnc|^^|mmol/l||><CR><ETX>7A<CR><LF>
<STX>6C|1|I|47C3 5000^0008 0000^Concentration > High System Limit|I<CR><ETX>7C<CR><LF>
<STX>7R|4|2028-9^LN^CO2 SerPl-sCnc|^^|mmol/l||<CR><ETX>E0<CR><LF>
<STX>0C|1|I|41CC 0000^0020 0080^False Rate (Unexpected Reaction Occurred)|I<CR><ETX>35<CR><LF>
<STX>10|2|||^* QUALITY CONTROL REPORT *: 4200D|20060401142819||||4||||4|||||||Q<CR><ETX>CD<CR><LF>
<STX>2C|1|I|^CHEMISTRY QC: 0|G<CR><ETX>2D<CR><LF>
<STX>3C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>10<CR><LF>
<STX>4R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>C8<CR><LF>
<STX>5R|2|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CA<CR><LF>
<STX>6R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>CE<CR><LF>
<STX>7R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|5|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D5<CR><LF>
<STX>1R|6|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>DB<CR><LF>
<STX>2R|7|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>E1<CR><LF>
<STX>3R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DB<CR><LF>
<STX>4R|9|^LEVEL 2: PRE|101|95 to 105|<CR><ETX>92<CR><LF>
<STX>5R|10|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>45<CR><LF>
<STX>6R|11|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|12|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>54<CR><LF>
<STX>0R|13|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>1R|14|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>52<CR><LF>
<STX>2R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>52<CR><LF>
<STX>3R|16|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>47<CR><LF>
<STX>4R|17|^LEVEL 2: 630 nm|101|95 to 105|<CR><ETX>6E<CR><LF>
<STX>5L|1|N<CR><ETX>05<CR><LF>
```


<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401145111<CR><ETX>E3<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^|Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1|^^^LIPID PANEL: 4600D|20060401144629||||7|||||7|||||||Q<CR><ETX>DC<CR><LF>
<STX>4C|1|I|^|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2093-3^^LN^Cholest SerPl-mCnc|ICT|mg/dl||<CR><ETX>CC<CR><LF>
<STX>2C|1|I|4353 0000^0200 0000^Icteric Limit Error|I<CR><ETX>56<CR><LF>
<STX>3R|2|2085-9^^LN^HDLc SerPl-mCnc|LIP|mg/dl||<CR><ETX>3C<CR><LF>
<STX>4C|1|I|4210 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>4C<CR><LF>
<STX>5R|3|2571-8^^LN^Trigl SerPl-mCnc|HEM|mg/dl||<CR><ETX>FA<CR><LF>
<STX>6C|1|I|42E6 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>71<CR><LF>
<STX>7R|4|9830-1^^LN^Cholest/HDLc SerPl-mRto|^^^||<CR><ETX>23<CR><LF>
<STX>0C|1|I|40BC CCCD^0000 0001^Chemistry Specific Error|I<CR><ETX>D6<CR><LF>
<STX>1R|5|13457-7^^LN^LDLc SerPl Calc-mCnc|152c|mg/dl||<CR><ETX>25<CR><LF>
<STX>2R|6|13458-5^^LN^VLDLc SerPl Calc-mCnc|23c|mg/dl||<CR><ETX>49<CR><LF>
<STX>30|2|^^^* QUALITY CONTROL REPORT *: 4600D|20060401144629||||7|||||7|||||||Q<CR><ETX>D2<CR><LF>
<STX>4C|1|I|^|^CHEMISTRY QC: 99|G<CR><ETX>49<CR><LF>
<STX>5C|2|I|^|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>12<CR><LF>
<STX>6R|1|^^^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CA<CR><LF>
<STX>7R|2|^^^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CC<CR><LF>
<STX>0R|3|^^^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D0<CR><LF>
<STX>1R|4|^^^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D3<CR><LF>
<STX>2R|5|^^^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D7<CR><LF>
<STX>3R|6|^^^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D5<CR><LF>
<STX>4R|7|^^^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DB<CR><LF>
<STX>5R|8|^^^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DD<CR><LF>
<STX>6R|9|^^^LEVEL 2: PRE|103|95 to 105|<CR><ETX>96<CR><LF>
<STX>7R|10|^^^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>47<CR><LF>
<STX>0R|11|^^^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4C<CR><LF>
<STX>1R|12|^^^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>56<CR><LF>
<STX>2R|13|^^^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>4C<CR><LF>
<STX>3R|14|^^^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4C<CR><LF>
<STX>4R|15|^^^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4C<CR><LF>
<STX>5R|16|^^^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>49<CR><LF>
<STX>6R|17|^^^LEVEL 2: 630 nm|103|95 to 105|<CR><ETX>72<CR><LF>
<STX>7L|1|N<CR><ETX>07<CR><LF>

*** Example (LIS to Abaxis) Query message ***

*** patient ID LEVEL IX, query results for LEVEL IX CONTROL rotor dumped below ***

```
<STX>1H|¥&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|LEVEL IX<CR><ETX>F5<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401150945<CR><ETX>EE<CR><LF>
<STX>2P|1|LEVEL IX||||U|||||^Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1||^METLYTE 8: 4300D||20060401150503||||0|||||0|||||||Q<CR><ETX>39<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2345-7^LN^Glucose SerPl-mCnc|HEM|mg/dl||<CR><ETX>CA<CR><LF>
<STX>2C|1|I|42BF 0000^OE00 0000^Hemolysis Limit Error|I<CR><ETX>82<CR><LF>
<STX>3R|2|3094-0^LN^BUN SerPl-mCnc|15|mg/dl||<CR><ETX>5F<CR><LF>
<STX>4R|3|2160-0^LN^Creat SerPl-mCnc|LIP|mg/dl||<CR><ETX>E3<CR><LF>
<STX>5C|1|I|3F66 6666^0600 0000^Lipemic Limit Error|I<CR><ETX>83<CR><LF>
<STX>6R|4|2157-6^LN^CK SerPl-cCnc|115|units/l||<CR><ETX>34<CR><LF>
<STX>7R|5|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F2<CR><LF>
<STX>0C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5C<CR><LF>
<STX>1R|6|2823-3^LN^Potassium SerPl-sCnc|4.0|mmol/l||<CR><ETX>FA<CR><LF>
<STX>2R|7|2075-0^LN^Chloride SerPl-sCnc|^^^|mmol/l||><CR><ETX>83<CR><LF>
<STX>3C|1|I|461C 4000^0008 0000^Concentration > High System Limit|I<CR><ETX>75<CR><LF>
<STX>4R|8|2028-9^LN^CO2 SerPl-sCnc|26|mmol/l||<CR><ETX>CF<CR><LF>
<STX>50|2|||^* QUALITY CONTROL REPORT *: 4300D||20060401150503||||0|||||0|||||||Q<CR><ETX>B7<CR><LF>
<STX>6C|1|I|^CHEMISTRY QC: 99|G<CR><ETX>4B<CR><LF>
<STX>7C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>14<CR><LF>
<STX>0R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CC<CR><LF>
<STX>1R|2|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CE<CR><LF>
<STX>2R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D2<CR><LF>
<STX>3R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>CD<CR><LF>
<STX>4R|5|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D1<CR><LF>
<STX>5R|6|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D7<CR><LF>
<STX>6R|7|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DD<CR><LF>
<STX>7R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DF<CR><LF>
<STX>0R|9|^LEVEL 2: PRE|100|95 to 105|<CR><ETX>95<CR><LF>
<STX>1R|10|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>49<CR><LF>
<STX>2R|11|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>3R|12|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>50<CR><LF>
<STX>4R|13|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>46<CR><LF>
<STX>5R|14|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>6R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4E<CR><LF>
<STX>7R|16|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>4B<CR><LF>
<STX>0R|17|^LEVEL 2: 630 nm|100|95 to 105|<CR><ETX>71<CR><LF>
<STX>1L|1|N<CR><ETX>09<CR><LF>
```

*** Example (LIS to Abaxis) Query message ***

*** patient ID ALL, query results for all rotors with results between 2006 04/01 14:00:19 & 2006 04/01 14:45:03 dumped below ***

```
<STX>1H|¥&|||Heisenberg|Union City|510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|ALL|||20060401140019|20060401144503<CR><ETX>78<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401143301<CR><ETX>E2<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^|Control|||||||<CR><ETX>21<CR><LF>
<STX>3O|1|||^ELECTROLYTE PANEL: 4200D|20060401142819||||4||||4|||||||Q<CR><ETX>AB<CR><LF>
<STX>4C|1|I|^INST QC: OK    CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>2C<CR><LF>
<STX>0C|5|I|^Control Renewal Required|G<CR><ETX>7D<CR><LF>
<STX>1R|1|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F0<CR><LF>
<STX>2C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5E<CR><LF>
<STX>3R|2|2823-3^LN^Potassium SerPl-sCnc|LIP|mmol/l||<CR><ETX>43<CR><LF>
<STX>4C|1|I|4096 6666^0600 0000^Lipemic Limit Error|I<CR><ETX>70<CR><LF>
<STX>5R|3|2075-0^LN^Chloride SerPl-sCnc|^^^|mmol/l||><CR><ETX>7A<CR><LF>
<STX>6C|1|I|47C3 5000^0008 0000^Concentration > High System Limit|I<CR><ETX>7C<CR><LF>
<STX>7R|4|2028-9^LN^CO2 SerPl-sCnc|^^^|mmol/l||<CR><ETX>E0<CR><LF>
<STX>0C|1|I|41CC 0000^0020 0080^False Rate (Unexpected Reaction Occurred)|I<CR><ETX>35<CR><LF>
<STX>1O|2|||^* QUALITY CONTROL REPORT *: 4200D|20060401142819||||4||||4|||||||Q<CR><ETX>CD<CR><LF>
<STX>2C|1|I|^CHEMISTRY QC: 0|G<CR><ETX>2D<CR><LF>
<STX>3C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>10<CR><LF>
<STX>4R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>C8<CR><LF>
<STX>5R|2|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CA<CR><LF>
<STX>6R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>CE<CR><LF>
<STX>7R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|5|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D5<CR><LF>
<STX>1R|6|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>DB<CR><LF>
<STX>2R|7|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>E1<CR><LF>
<STX>3R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DB<CR><LF>
<STX>4R|9|^LEVEL 2: PRE|101|95 to 105|<CR><ETX>92<CR><LF>
<STX>5R|10|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>45<CR><LF>
<STX>6R|11|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|12|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>54<CR><LF>
<STX>0R|13|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>1R|14|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>52<CR><LF>
<STX>2R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>52<CR><LF>
<STX>3R|16|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>47<CR><LF>
<STX>4R|17|^LEVEL 2: 630 nm|101|95 to 105|<CR><ETX>6E<CR><LF>
<STX>5L|1|N<CR><ETX>05<CR><LF>
```

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401144000<CR><ETX>DF<CR><LF>
<STX>2P|1|5|44444|33333|||19530225|M|||789-012-3456||88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^GENERAL CHEMISTRY 12: 1700D||20060401143518|||5|||||5|||||||Q<CR><ETX>29<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: 24 LIP: 166 ICT: 1|G<CR><ETX>17<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|1751-7^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>
<STX>2C|1|I|408C CCCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>
<STX>3R|2|6768-6^LN^ALP SerPl-cCnc|55|units/l|53 to 128|N<CR><ETX>D1<CR><LF>
<STX>4R|3|1742-6^LN^ALT SerPl-cCnc|HEM|units/l|10 to 47|<CR><ETX>B5<CR><LF>
<STX>5C|1|I|41E4 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>6D<CR><LF>
<STX>6R|4|1798-8^LN^Amylase SerPl-cCnc|56|units/l|14 to 97|N<CR><ETX>98<CR><LF>
<STX>7R|5|1920-8^LN^AST SerPl-cCnc|LIP|units/l|11 to 38|<CR><ETX>CD<CR><LF>
<STX>0C|1|I|41C4 0000^0600 0100^Lipemic Limit Error|I<CR><ETX>66<CR><LF>
<STX>1R|6|1975-2^LN^Bilirub SerPl-mCnc|0.9|mg/dl|0.2 to 1.6|N<CR><ETX>1C<CR><LF>
<STX>2R|7|3094-0^LN^BUN SerPl-mCnc|15|mg/dl|7 to 22|N<CR><ETX>77<CR><LF>
<STX>3R|8|17861-6^LN^Calcium SerPl-mCnc|ICT|mg/dl|8.0 to 10.3|<CR><ETX>70<CR><LF>
<STX>4C|1|I|4112 6666^0200 0002^Icteric Limit Error|I<CR><ETX>63<CR><LF>
<STX>5R|9|2093-3^LN^Cholest SerPl-mCnc|150|mg/dl|100 to 200|N<CR><ETX>1A<CR><LF>
<STX>6R|10|2160-0^LN^Creat SerPl-mCnc|0.9|mg/dl|0.6 to 1.2|N<CR><ETX>5B<CR><LF>
<STX>7R|11|2345-7^LN^Glucose SerPl-mCnc|96|mg/dl|73 to 118|N<CR><ETX>03<CR><LF>
<STX>0R|12|2885-2^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>C3<CR><LF>
<STX>10|2|||^* QUALITY CONTROL REPORT *: 1700D||20060401143518|||5|||||5|||||||Q<CR><ETX>CE<CR><LF>
<STX>2C|1|I|^CHEMISTRY QC: 0|G<CR><ETX>2D<CR><LF>
<STX>3C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>10<CR><LF>
<STX>4R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>C8<CR><LF>
<STX>5R|2|^LEVEL 1: IQC 2|102|90 to 110|<CR><ETX>CB<CR><LF>
<STX>6R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>CE<CR><LF>
<STX>7R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|5|^LEVEL 1: IQC 5|101|90 to 110|<CR><ETX>D3<CR><LF>
<STX>1R|6|^LEVEL 1: IQC 6|101|90 to 110|<CR><ETX>D6<CR><LF>
<STX>2R|7|^LEVEL 1: IQC 7|106|90 to 110|<CR><ETX>DE<CR><LF>
<STX>3R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DB<CR><LF>
<STX>4R|9|^LEVEL 2: PRE|102|95 to 105|<CR><ETX>93<CR><LF>
<STX>5R|10|^LEVEL 2: 340 nm|99|95 to 105|<CR><ETX>46<CR><LF>
<STX>6R|11|^LEVEL 2: 405 nm|100|95 to 105|<CR><ETX>69<CR><LF>
<STX>7R|12|^LEVEL 2: 467 nm|100|95 to 105|<CR><ETX>73<CR><LF>
<STX>0R|13|^LEVEL 2: 500 nm|100|95 to 105|<CR><ETX>69<CR><LF>
<STX>1R|14|^LEVEL 2: 515 nm|100|95 to 105|<CR><ETX>71<CR><LF>
<STX>2R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>52<CR><LF>
<STX>3R|16|^LEVEL 2: 600 nm|98|95 to 105|<CR><ETX>48<CR><LF>
<STX>4R|17|^LEVEL 2: 630 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>5L|1|N<CR><ETX>05<CR><LF>

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401144524<CR><ETX>EA<CR><LF>
<STX>2P|1|6|44444|33333|||19530225|F|||789-012-3456|88888^53 Yrs.^|Patient|||||||66666<CR><ETX>21<CR><LF>
<STX>30|1|||^^^LIVER PANEL PLUS: 1100D|20060401144042||||6|||||||Q<CR><ETX>3B<CR><LF>
<STX>4C|1|I|^|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|1751-7^^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>
<STX>2C|1|I|408C CCCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>
<STX>3R|2|6768-6^^LN^ALP SerPl-cCnc|55|units/l|42 to 141|N<CR><ETX>CA<CR><LF>
<STX>4R|3|1742-6^^LN^ALT SerPl-cCnc|LIP|units/l|10 to 47|<CR><ETX>C0<CR><LF>
<STX>5C|1|I|41E4 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>64<CR><LF>
<STX>6R|4|1798-8^^LN^Amylase SerPl-cCnc|56|units/l|14 to 97|N<CR><ETX>98<CR><LF>
<STX>7R|5|1920-8^^LN^AST SerPl-cCnc|<5*|units/l|11 to 38|<CR><ETX>BF<CR><LF>
<STX>0C|1|I|41C4 0000^0001 0000^Concentration < Low Dynamic Limit|I<CR><ETX>5D<CR><LF>
<STX>1R|6|1975-2^^LN^Bilirub SerPl-mCnc|0.9|mg/dl|0.2 to 1.6|N<CR><ETX>1C<CR><LF>
<STX>2R|7|2324-2^^LN^GGT SerPl-cCnc|35|units/l|5 to 65|N<CR><ETX>69<CR><LF>
<STX>3R|8|2885-2^^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>93<CR><LF>
<STX>40|2|||^|^* QUALITY CONTROL REPORT *: 1100D|20060401144042||||6|||||6|||||||Q<CR><ETX>BE<CR><LF>
<STX>5C|1|I|^|^CHEMISTRY QC: 0|G<CR><ETX>28<CR><LF>
<STX>6C|2|I|^|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>13<CR><LF>
<STX>7R|1|^|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CB<CR><LF>
<STX>0R|2|^|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CD<CR><LF>
<STX>1R|3|^|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>2R|4|^|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D4<CR><LF>
<STX>3R|5|^|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D0<CR><LF>
<STX>4R|6|^|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D6<CR><LF>
<STX>5R|7|^|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DC<CR><LF>
<STX>6R|8|^|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DE<CR><LF>
<STX>7R|9|^|^LEVEL 2: PRE|103|95 to 105|<CR><ETX>97<CR><LF>
<STX>0R|10|^|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>48<CR><LF>
<STX>1R|11|^|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4D<CR><LF>
<STX>2R|12|^|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>57<CR><LF>
<STX>3R|13|^|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>45<CR><LF>
<STX>4R|14|^|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4D<CR><LF>
<STX>5R|15|^|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4D<CR><LF>
<STX>6R|16|^|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|17|^|^LEVEL 2: 630 nm|103|95 to 105|<CR><ETX>73<CR><LF>
<STX>0L|1|N<CR><ETX>08<CR><LF>

*** Example (LIS to Abaxis) Query message ***

*** patient ID ALL, query results for all rotors with results ending 2006 04/01 14:28:19 dumped below ***

```
<STX>1H|¥&|||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|ALL|||20060401142819<CR><ETX>BF<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401143301<CR><ETX>E2<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^|Control|||||||<CR><ETX>21<CR><LF>
<STX>3O|1|||^ELECTROLYTE PANEL: 4200D|20060401142819||||4||||4|||||||Q<CR><ETX>AB<CR><LF>
<STX>4C|1|I|^INST QC: OK    CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>2C<CR><LF>
<STX>0C|5|I|^Control Renewal Required|G<CR><ETX>7D<CR><LF>
<STX>1R|1|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F0<CR><LF>
<STX>2C|1|I|4308 8000^0200 0000^Icteric Limit Error|I<CR><ETX>5E<CR><LF>
<STX>3R|2|2823-3^LN^Potassium SerPl-sCnc|LIP|mmol/l||<CR><ETX>43<CR><LF>
<STX>4C|1|I|4096 6666^0600 0000^Lipemic Limit Error|I<CR><ETX>70<CR><LF>
<STX>5R|3|2075-0^LN^Chloride SerPl-sCnc|^^|mmol/l||><CR><ETX>7A<CR><LF>
<STX>6C|1|I|47C3 5000^0008 0000^Concentration > High System Limit|I<CR><ETX>7C<CR><LF>
<STX>7R|4|2028-9^LN^CO2 SerPl-sCnc|^^|mmol/l||<CR><ETX>E0<CR><LF>
<STX>0C|1|I|41CC 0000^0020 0080^False Rate (Unexpected Reaction Occurred)|I<CR><ETX>35<CR><LF>
<STX>10|2|||^* QUALITY CONTROL REPORT *: 4200D|20060401142819||||4||||4|||||||Q<CR><ETX>CD<CR><LF>
<STX>2C|1|I|^CHEMISTRY QC: 0|G<CR><ETX>2D<CR><LF>
<STX>3C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>10<CR><LF>
<STX>4R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>C8<CR><LF>
<STX>5R|2|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CA<CR><LF>
<STX>6R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>CE<CR><LF>
<STX>7R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D1<CR><LF>
<STX>0R|5|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D5<CR><LF>
<STX>1R|6|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>DB<CR><LF>
<STX>2R|7|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>E1<CR><LF>
<STX>3R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DB<CR><LF>
<STX>4R|9|^LEVEL 2: PRE|101|95 to 105|<CR><ETX>92<CR><LF>
<STX>5R|10|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>45<CR><LF>
<STX>6R|11|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>7R|12|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>54<CR><LF>
<STX>0R|13|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>4A<CR><LF>
<STX>1R|14|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>52<CR><LF>
<STX>2R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>52<CR><LF>
<STX>3R|16|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>47<CR><LF>
<STX>4R|17|^LEVEL 2: 630 nm|101|95 to 105|<CR><ETX>6E<CR><LF>
<STX>5L|1|N<CR><ETX>05<CR><LF>
```

*** Example (LIS to Abaxis) Query message ***

*** patient ID ALL, query results for all rotors with results beginning 2006 04/01 15:05:03 dumped below ***

```
<STX>1H|¥&||Heisenberg|Union City|510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|ALL|||20060401150503<CR><ETX>3E<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401150945<CR><ETX>EE<CR><LF>
<STX>2P|1|LEVEL IX||||U|||||^Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1||^METLYTE 8: 4300D|20060401150503|||0|||||0|||||||Q<CR><ETX>39<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^G<CR><ETX>19<CR><LF>
<STX>1R|1|2345-7^LN^Glucose SerPl-mCnc|HEM|mg/dl||<CR><ETX>CA<CR><LF>
<STX>2C|1|I|42BF 0000^OE00 0000^Hemolysis Limit Error|I<CR><ETX>82<CR><LF>
<STX>3R|2|3094-0^LN^BUN SerPl-mCnc|15|mg/dl||<CR><ETX>5F<CR><LF>
<STX>4R|3|2160-0^LN^Creat SerPl-mCnc|LIP|mg/dl||<CR><ETX>E3<CR><LF>
<STX>5C|1|I|3F66 6666^O600 0000^Lipemic Limit Error|I<CR><ETX>83<CR><LF>
<STX>6R|4|2157-6^LN^CK SerPl-cCnc|115|units/l||<CR><ETX>34<CR><LF>
<STX>7R|5|2951-2^LN^Sodium SerPl-sCnc|ICT|mmol/l||<CR><ETX>F2<CR><LF>
<STX>0C|1|I|4308 8000^O200 0000^Icteric Limit Error|I<CR><ETX>5C<CR><LF>
<STX>1R|6|2823-3^LN^Potassium SerPl-sCnc|4.0|mmol/l||<CR><ETX>FA<CR><LF>
<STX>2R|7|2075-0^LN^Chloride SerPl-sCnc|~~~|mmol/l||><CR><ETX>83<CR><LF>
<STX>3C|1|I|461C 4000^O008 0000^Concentration > High System Limit|I<CR><ETX>75<CR><LF>
<STX>4R|8|2028-9^LN^CO2 SerPl-sCnc|26|mmol/l||<CR><ETX>CF<CR><LF>
<STX>50|2|||^* QUALITY CONTROL REPORT *: 4300D|20060401150503|||0|||||0|||||||Q<CR><ETX>B7<CR><LF>
<STX>6C|1|I|^CHEMISTRY QC: 99|G<CR><ETX>4B<CR><LF>
<STX>7C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>14<CR><LF>
<STX>0R|1|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CC<CR><LF>
<STX>1R|2|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CE<CR><LF>
<STX>2R|3|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D2<CR><LF>
<STX>3R|4|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>CD<CR><LF>
<STX>4R|5|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D1<CR><LF>
<STX>5R|6|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D7<CR><LF>
<STX>6R|7|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DD<CR><LF>
<STX>7R|8|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DF<CR><LF>
<STX>0R|9|^LEVEL 2: PRE|100|95 to 105|<CR><ETX>95<CR><LF>
<STX>1R|10|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>49<CR><LF>
<STX>2R|11|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>3R|12|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>50<CR><LF>
<STX>4R|13|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>46<CR><LF>
<STX>5R|14|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4E<CR><LF>
<STX>6R|15|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4E<CR><LF>
<STX>7R|16|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>4B<CR><LF>
<STX>0R|17|^LEVEL 2: 630 nm|100|95 to 105|<CR><ETX>71<CR><LF>
```

<STX>1L|1|N<CR><ETX>09<CR><LF>

*** Example (LIS to Abaxis) Query message ***

*** patient ID ALL, query results for all rotors with results between 2006 04/01 14:41:20 & 2006 04/01 14:58:02 dumped below ***

<STX>1H|ÿ^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2Q|1|ALL|||20060401144120|20060401145802<CR><ETX>6F<CR><LF>

<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|ÿ^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401145111<CR><ETX>E3<CR><LF>

<STX>2P|1|LEVEL V|||||U|||||^|^|Control|||||||<CR><ETX>21<CR><LF>

<STX>30|1|||^|^|^LIPID PANEL: 4600D|20060401144629||||7|||||7|||||||Q<CR><ETX>DC<CR><LF>

<STX>4C|1|I|^|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>

<STX>5C|2|I|^|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>

<STX>6C|3|I|^|^|G<CR><ETX>1D<CR><LF>

<STX>7C|4|I|^|^|G<CR><ETX>1F<CR><LF>

<STX>0C|5|I|^|^|G<CR><ETX>19<CR><LF>

<STX>1R|1|2093-3^LN^Cholest SerPl-mCnc|ICT|mg/dl||<CR><ETX>CC<CR><LF>

<STX>2C|1|I|4353 0000^0200 0000^Icteric Limit Error|I<CR><ETX>56<CR><LF>

<STX>3R|2|2085-9^LN^HDLc SerPl-mCnc|LIP|mg/dl||<CR><ETX>3C<CR><LF>

<STX>4C|1|I|4210 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>4C<CR><LF>

<STX>5R|3|2571-8^LN^Trigl SerPl-mCnc|HEM|mg/dl||<CR><ETX>FA<CR><LF>

<STX>6C|1|I|42E6 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>71<CR><LF>

<STX>7R|4|9830-1^LN^Cholest/HDLc SerPl-mRto|^|^|^||<CR><ETX>23<CR><LF>

<STX>0C|1|I|40BC CCCD^0000 0001^Chemistry Specific Error|I<CR><ETX>D6<CR><LF>

<STX>1R|5|13457-7^LN^LDLc SerPl Calc-mCnc|152c|mg/dl||<CR><ETX>25<CR><LF>

<STX>2R|6|13458-5^LN^VLDLc SerPl Calc-mCnc|23c|mg/dl||<CR><ETX>49<CR><LF>

<STX>30|2|||^|^* QUALITY CONTROL REPORT *: 4600D|20060401144629||||7|||||7|||||||Q<CR><ETX>D2<CR><LF>

<STX>4C|1|I|^|^CHEMISTRY QC: 99|G<CR><ETX>49<CR><LF>

<STX>5C|2|I|^|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>12<CR><LF>

<STX>6R|1|^|^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CA<CR><LF>

<STX>7R|2|^|^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CC<CR><LF>

<STX>0R|3|^|^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D0<CR><LF>

<STX>1R|4|^|^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D3<CR><LF>

<STX>2R|5|^|^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D7<CR><LF>

<STX>3R|6|^|^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D5<CR><LF>

<STX>4R|7|^|^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DB<CR><LF>

<STX>5R|8|^|^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DD<CR><LF>

<STX>6R|9|^|^LEVEL 2: PRE|103|95 to 105|<CR><ETX>96<CR><LF>

<STX>7R|10|^|^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>47<CR><LF>

<STX>0R|11|^|^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4C<CR><LF>

<STX>1R|12|^|^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>56<CR><LF>

<STX>2R|13|^|^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>4C<CR><LF>

<STX>3R|14|^|^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4C<CR><LF>

<STX>4R|15|^|^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4C<CR><LF>

<STX>5R|16|^|^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>49<CR><LF>

<STX>6R|17|^|^LEVEL 2: 630 nm|103|95 to 105|<CR><ETX>72<CR><LF>

<STX>7L|1|N<CR><ETX>07<CR><LF>

<STX>1H|¥&|||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D||||||D|E 1394-97|20060401145704<CR><ETX>EB<CR><LF>

<STX>2P|1|8|44444|33333|||19530225|U|||789-012-3456|88888^53 Yrs.^Patient|||||||66666<CR><ETX>21<CR><LF>

<STX>30|1|||^^^HEPATIC FUNCTION PANEL: 4700D|20060401145222||||8|||||8|||||||Q<CR><ETX>E7<CR><LF>

<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>

<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>

<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>

<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>

<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>

<STX>1R|1|1751-7^LN^Alb SerPl-mCnc|ICT|g/dl|3.3 to 5.5|<CR><ETX>EF<CR><LF>

<STX>2C|1|I|408C CCCD^0200 0000^Icteric Limit Error|I<CR><ETX>B3<CR><LF>

<STX>3R|2|6768-6^LN^ALP SerPl-cCnc|55|units/l|42 to 141|N<CR><ETX>CA<CR><LF>

<STX>4R|3|1742-6^LN^ALT SerPl-cCnc|ICT|units/l|10 to 47|<CR><ETX>BB<CR><LF>

<STX>5C|1|I|41E8 0000^0200 0000^Icteric Limit Error|I<CR><ETX>64<CR><LF>

<STX>6R|4|1920-8^LN^AST SerPl-cCnc|^^^|units/l|11 to 38|<CR><ETX>60<CR><LF>

<STX>7C|1|I|41C8 0000^0000 0108^Chemistry Specific Error|I<CR><ETX>87<CR><LF>

<STX>0R|5|1968-7^LN^Bilirub Direct SerPl-mCnc|0.2|mg/dl|0.0 to 0.3|N<CR><ETX>8F<CR><LF>

<STX>1R|6|1975-2^LN^Bilirub SerPl-mCnc|^^^|mg/dl|0.2 to 1.6|<<CR><ETX>ED<CR><LF>

<STX>2C|1|I|3F66 6666^0004 0000^Concentration < Low System Limit|I<CR><ETX>43<CR><LF>

<STX>3R|7|2885-2^LN^Prot SerPl-mCnc|7.3|g/dl|6.4 to 8.1|N<CR><ETX>92<CR><LF>

<STX>40|2|||^* QUALITY CONTROL REPORT *: 4700D|20060401145222||||8|||||8|||||||Q<CR><ETX>CC<CR><LF>

<STX>5C|1|I|^CHEMISTRY QC: 99|G<CR><ETX>4A<CR><LF>

<STX>6C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>13<CR><LF>

<STX>7R|1|^LEVEL 1: IQC 1|102||90 to 110|<CR><ETX>CB<CR><LF>

<STX>0R|2|^LEVEL 1: IQC 2|101||90 to 110|<CR><ETX>CD<CR><LF>

<STX>1R|3|^LEVEL 1: IQC 3|102||90 to 110|<CR><ETX>D1<CR><LF>

<STX>2R|4|^LEVEL 1: IQC 4|102||90 to 110|<CR><ETX>D4<CR><LF>

<STX>3R|5|^LEVEL 1: IQC 5|103||90 to 110|<CR><ETX>D0<CR><LF>

<STX>4R|6|^LEVEL 1: IQC 6|106||90 to 110|<CR><ETX>D6<CR><LF>

<STX>5R|7|^LEVEL 1: IQC 7|109||90 to 110|<CR><ETX>DC<CR><LF>

<STX>6R|8|^LEVEL 1: IQC 8|108||90 to 110|<CR><ETX>DE<CR><LF>

<STX>7R|9|^LEVEL 2: PRE|103||95 to 105|<CR><ETX>97<CR><LF>

<STX>0R|10|^LEVEL 2: 340 nm|98||95 to 105|<CR><ETX>48<CR><LF>

<STX>1R|11|^LEVEL 2: 405 nm|99||95 to 105|<CR><ETX>4D<CR><LF>

<STX>2R|12|^LEVEL 2: 467 nm|99||95 to 105|<CR><ETX>57<CR><LF>

<STX>3R|13|^LEVEL 2: 500 nm|99||95 to 105|<CR><ETX>45<CR><LF>

<STX>4R|14|^LEVEL 2: 515 nm|99||95 to 105|<CR><ETX>4D<CR><LF>

<STX>5R|15|^LEVEL 2: 550 nm|98||95 to 105|<CR><ETX>4D<CR><LF>

<STX>6R|16|^LEVEL 2: 600 nm|97||95 to 105|<CR><ETX>4A<CR><LF>

<STX>7R|17|^LEVEL 2: 630 nm|103||95 to 105|<CR><ETX>73<CR><LF>

<STX>0L|1|N<CR><ETX>08<CR><LF>

*** Example (LIS to Abaxis) Query message ***

*** patient ID LEVEL 5, query results for rotors with results between 2006 04/01 14:41:20 & 2006 04/01 14:58:02 dumped below ***

```
<STX>1H|¥&||Heisenberg|Union City|510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|LEVEL V|||20060401144120|20060401145802<CR><ETX>D1<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401145111<CR><ETX>E3<CR><LF>
<STX>2P|1|LEVEL V|||||U|||||^Control|||||||<CR><ETX>21<CR><LF>
<STX>30|1|^^^LIPID PANEL: 4600D|20060401144629|||7|||||7|||||||Q<CR><ETX>DC<CR><LF>
<STX>4C|1|I|^INST QC: OK CHEM QC: OK|G<CR><ETX>1F<CR><LF>
<STX>5C|2|I|^HEM: -32 LIP: 210 ICT: 5|G<CR><ETX>1D<CR><LF>
<STX>6C|3|I|^|G<CR><ETX>1D<CR><LF>
<STX>7C|4|I|^|G<CR><ETX>1F<CR><LF>
<STX>0C|5|I|^|G<CR><ETX>19<CR><LF>
<STX>1R|1|2093-3^LN^Cholest SerPl-mCnc|ICT|mg/dl||<CR><ETX>CC<CR><LF>
<STX>2C|1|I|4353 0000^0200 0000^Icteric Limit Error|I<CR><ETX>56<CR><LF>
<STX>3R|2|2085-9^LN^HDLc SerPl-mCnc|LIP|mg/dl||<CR><ETX>3C<CR><LF>
<STX>4C|1|I|4210 0000^0600 0000^Lipemic Limit Error|I<CR><ETX>4C<CR><LF>
<STX>5R|3|2571-8^LN^Trigl SerPl-mCnc|HEM|mg/dl||<CR><ETX>FA<CR><LF>
<STX>6C|1|I|42E6 0000^0E00 0000^Hemolysis Limit Error|I<CR><ETX>71<CR><LF>
<STX>7R|4|9830-1^LN^Cholest/HDLc SerPl-mRto|^^^||<CR><ETX>23<CR><LF>
<STX>0C|1|I|40BC CCCD^0000 0001^Chemistry Specific Error|I<CR><ETX>D6<CR><LF>
<STX>1R|5|13457-7^LN^LDLc SerPl Calc-mCnc|152c|mg/dl||<CR><ETX>25<CR><LF>
<STX>2R|6|13458-5^LN^VLDLc SerPl Calc-mCnc|23c|mg/dl||<CR><ETX>49<CR><LF>
<STX>30|2|^^^* QUALITY CONTROL REPORT *: 4600D|20060401144629|||7|||||7|||||||Q<CR><ETX>D2<CR><LF>
<STX>4C|1|I|^CHEMISTRY QC: 99|G<CR><ETX>49<CR><LF>
<STX>5C|2|I|^ACCEPTABLE MINIMUM: 50|G<CR><ETX>12<CR><LF>
<STX>6R|1|^^^LEVEL 1: IQC 1|102|90 to 110|<CR><ETX>CA<CR><LF>
<STX>7R|2|^^^LEVEL 1: IQC 2|101|90 to 110|<CR><ETX>CC<CR><LF>
<STX>0R|3|^^^LEVEL 1: IQC 3|102|90 to 110|<CR><ETX>D0<CR><LF>
<STX>1R|4|^^^LEVEL 1: IQC 4|102|90 to 110|<CR><ETX>D3<CR><LF>
<STX>2R|5|^^^LEVEL 1: IQC 5|103|90 to 110|<CR><ETX>D7<CR><LF>
<STX>3R|6|^^^LEVEL 1: IQC 6|106|90 to 110|<CR><ETX>D5<CR><LF>
<STX>4R|7|^^^LEVEL 1: IQC 7|109|90 to 110|<CR><ETX>DB<CR><LF>
<STX>5R|8|^^^LEVEL 1: IQC 8|108|90 to 110|<CR><ETX>DD<CR><LF>
<STX>6R|9|^^^LEVEL 2: PRE|103|95 to 105|<CR><ETX>96<CR><LF>
<STX>7R|10|^^^LEVEL 2: 340 nm|98|95 to 105|<CR><ETX>47<CR><LF>
<STX>0R|11|^^^LEVEL 2: 405 nm|99|95 to 105|<CR><ETX>4C<CR><LF>
<STX>1R|12|^^^LEVEL 2: 467 nm|99|95 to 105|<CR><ETX>56<CR><LF>
<STX>2R|13|^^^LEVEL 2: 500 nm|99|95 to 105|<CR><ETX>4C<CR><LF>
<STX>3R|14|^^^LEVEL 2: 515 nm|99|95 to 105|<CR><ETX>4C<CR><LF>
<STX>4R|15|^^^LEVEL 2: 550 nm|98|95 to 105|<CR><ETX>4C<CR><LF>
<STX>5R|16|^^^LEVEL 2: 600 nm|97|95 to 105|<CR><ETX>49<CR><LF>
<STX>6R|17|^^^LEVEL 2: 630 nm|103|95 to 105|<CR><ETX>72<CR><LF>
<STX>7L|1|N<CR><ETX>07<CR><LF>
```

***** Example (LIS to Abaxis) Query message *****

***** patient ID NULL ('¥0'), no query results message sent below for ID [NULL ('¥0')] for rotors with results between 2006 04/01 14:00:00 & 2006 04/01 14:27:28 *****

```
<STX>1H|¥^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|||||20060401140000|20060401142728<CR><ETX>B0<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155140<CR><ETX>E6<CR><LF>
<STX>2P|1|<CR><ETX>BB<CR><LF>
<STX>30|1|||||||||||||||||||||Z<CR><ETX>39<CR><LF>
<STX>4C|1|I|^:20060401140000:20060401142728:No Results for this Query|G<CR><ETX>6E<CR><LF>
<STX>5L|1|I<CR><ETX>03<CR><LF>
```

***** Example (LIS to Abaxis) Query message *****

***** patient ID NULL ('¥0'), no query results message sent below for ID [NULL ('¥0')] *****

```
<STX>1H|¥^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1<CR><ETX>40<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155200<CR><ETX>E3<CR><LF>
<STX>2P|1|<CR><ETX>BB<CR><LF>
<STX>30|1|||||||||||||||||||||Z<CR><ETX>39<CR><LF>
<STX>4C|1|I|^::No Results for this Query|G<CR><ETX>E6<CR><LF>
<STX>5L|1|I<CR><ETX>03<CR><LF>
```

***** Example (LIS to Abaxis) Query message *****

***** patient ID 11111, no query results message sent below for ID 11111 *****

```
<STX>1H|¥^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>
<STX>2Q|1|11111<CR><ETX>B1<CR><LF>
<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155221<CR><ETX>E6<CR><LF>
<STX>2P|1|11111<CR><ETX>B0<CR><LF>
<STX>30|1|||||||||||||||||||||Z<CR><ETX>39<CR><LF>
<STX>4C|1|I|^11111::No Results for this Query|G<CR><ETX>DB<CR><LF>
<STX>5L|1|I<CR><ETX>03<CR><LF>
```

*** patient ID LEVEL 5, no query results message sent below for ID LEVEL 5 ***

<STX>20|1|LEVEL 5<CR><ETX>89<CR><LF>

<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>2P|1|LEVEL 5<CR><ETX>88<CR><LF>

```
<STX>30 1|||||||Z<CR><ETX>39<CR><LF>
```

```
<STX>4C|1|I|^LEVEL 5:::No Results for this Query|G<CR><ETX>B3<CR><LF>
```

<STX>5L|1|I<CR><ETX>03<CR><LF>

*** patient ID level 9, no query results message sent below for ID level 9 ***

<STX>2Q|1|level 9<CR><ETX>2D<CR><LF>

<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>2P|1|LEVEL 9<CR><ETX>8C<CR><LF>

```
<STX>30|1|||Z<CR><ETX>39<CR><LF>
```

```
<STX>4C|1|I|^LEVEL 9::No Results for this Query|G<CR><ETX>B7<CR><LF>
```

<STX>5L|1|I<CR><ETX>03<CR><LF>

*** patient ID 999999999999999999, ^Q^Received Field Size Too Big & no query results messages sent below for ID 999999999999999999 ***

<STX>2Q|1|9999999999999999999<CR><ETX>30<CR><LF>

<STX>3L|1<CR><ETX>3C<CR><LF>

```
<STX>2C|1|I|^QReceived Field Size Too Big|G<CR><ETX>AD<CR><LF>
```

<STX>3L|1|N<CR><ETX>06<CR><LF>

```
<STX>2P|1|9999999999999999<CR><ETX>D9<CR><LF>
```

```
<STX>30|1|||Z<CR><ETX>39<CR><LF>
```

```
<STX>4C|1|I|^9999999999999999::No Results for this Query|G<CR><ETX>04<CR><LF>
```

<STX>5L|1|I<CR><ETX>03<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID 12345678901234, ^Q^Received BAD Q Record Sequence # & no query results messages sent below for ID 12345678901234 *****

<STX>1H|¥^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2Q|2|12345678901234<CR><ETX>94<CR><LF>

<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155514<CR><ETX>EB<CR><LF>

<STX>2C|1|I|^Q^Received BAD Q Record Sequence #|G<CR><ETX>DD<CR><LF>

<STX>3L|1|N<CR><ETX>06<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155516<CR><ETX>ED<CR><LF>

<STX>2P|1|12345678901234<CR><ETX>92<CR><LF>

<STX>30|1|||||||||||||||||||||Z<CR><ETX>39<CR><LF>

<STX>4C|1|I|^12345678901234::No Results for this Query|G<CR><ETX>BD<CR><LF>

<STX>5L|1|I<CR><ETX>03<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID 12345678901234, ^L^Received BAD L Record Sequence # & no query results messages sent below for ID 12345678901234 *****

<STX>1H|¥^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2Q|1|12345678901234<CR><ETX>93<CR><LF>

<STX>3L|2<CR><ETX>3D<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155554<CR><ETX>EF<CR><LF>

<STX>2C|1|I|^L^Received BAD L Record Sequence #|G<CR><ETX>D3<CR><LF>

<STX>3L|1|N<CR><ETX>06<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155556<CR><ETX>F1<CR><LF>

<STX>2P|1|12345678901234<CR><ETX>92<CR><LF>

<STX>30|1|||||||||||||||||||||Z<CR><ETX>39<CR><LF>

<STX>4C|1|I|^12345678901234::No Results for this Query|G<CR><ETX>BD<CR><LF>

<STX>5L|1|I<CR><ETX>03<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID 12345678901234, no query results message sent below for ID 12345678901234 *****

<STX>1H|¥^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2Q|1|12345678901234<CR><ETX>93<CR><LF>

<STX>3L|1|x<CR><ETX>30<CR><LF>

<STX>1H|¥^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155642<CR><ETX>ED<CR><LF>

<STX>2P|1|12345678901234<CR><ETX>92<CR><LF>

<STX>30|1|||||||||||||||||||||Z<CR><ETX>39<CR><LF>

<STX>4C|1|I|^12345678901234::No Results for this Query|G<CR><ETX>BD<CR><LF>

<STX>5L|1|I<CR><ETX>03<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID 12345678901234, ^S^Received Record Out of Sequence & no query results messages sent below for ID 12345678901234 *****

<STX>1Q|1|12345678901234<CR><ETX>92<CR><LF>

<STX>2H|^&||Heisenberg|Union City||510-675-6500<CR><ETX>9F<CR><LF>

<STX>3Q|1|12345678901234<CR><ETX>94<CR><LF>

<STX>4L|1|x<CR><ETX>31<CR><LF>

<STX>5S<CR><ETX>98<CR><LF>

<STX>1H|^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155715<CR><ETX>EE<CR><LF>

<STX>2C|1|I|^S^Received Record Out of Sequence|G<CR><ETX>91<CR><LF>

<STX>3L|1|N<CR><ETX>06<CR><LF>

<STX>1H|^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155716<CR><ETX>EF<CR><LF>

<STX>2P|1|12345678901234<CR><ETX>92<CR><LF>

<STX>30|1|||||||||||||||||Z<CR><ETX>39<CR><LF>

<STX>4C|1|I|^12345678901234:::No Results for this Query|G<CR><ETX>BD<CR><LF>

<STX>5L|1|I<CR><ETX>03<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** ^^Received Message NOT Supported message sent below *****

<STX>1H|^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2Q|1|12345678901234<CR><ETX>93<CR><LF>

<STX>3R<CR><ETX>95<CR><LF>

<STX>4L|1<CR><ETX>3D<CR><LF>

<STX>1H|^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401155930<CR><ETX>ED<CR><LF>

<STX>2C|1|I|^Received Message NOT Supported|G<CR><ETX>F5<CR><LF>

<STX>3L|1|N<CR><ETX>06<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** ^X^Received Invalid Record Type message sent below *****

<STX>1H|^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2X<CR><ETX>9A<CR><LF>

<STX>3L|1<CR><ETX>3C<CR><LF>

<STX>1H|^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401160011<CR><ETX>DF<CR><LF>

<STX>2C|1|I|^X^Received Invalid Record Type|G<CR><ETX>99<CR><LF>

<STX>3L|1|N<CR><ETX>06<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID 12345678901234, ^X^Received Invalid Record Type message sent below *****

<STX>1H|^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2Q|1|12345678901234<CR><ETX>93<CR><LF>

<STX>3X<CR><ETX>9B<CR><LF>

<STX>4L|1<CR><ETX>3D<CR><LF>

<STX>1H|^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401160033<CR><ETX>E3<CR><LF>

<STX>2C|1|I|^X^Received Invalid Record Type|G<CR><ETX>99<CR><LF>

<STX>3L|1|N<CR><ETX>06<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID 12345678901234, ^X^Received Invalid Record Type message sent below *****

<STX>1H|^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2X<CR><ETX>9A<CR><LF>

<STX>3Q|1|12345678901234<CR><ETX>94<CR><LF>

<STX>4L|1<CR><ETX>3D<CR><LF>

<STX>1H|^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401160047<CR><ETX>E8<CR><LF>

<STX>2C|1|I|^X^Received Invalid Record Type|G<CR><ETX>99<CR><LF>

<STX>3L|1|N<CR><ETX>06<CR><LF>

***** Example (LIS to Abaxis) Query message *****

***** patient ID 12345678901234, ^^Received Message NOT Supported message sent below *****

<STX>1H|^&||Heisenberg|Union City||510-675-6500<CR><ETX>9E<CR><LF>

<STX>2Q|1|12345678901234<CR><ETX>93<CR><LF>

<STX>3Q|1|12345678901234<CR><ETX>94<CR><LF>

<STX>4L|1<CR><ETX>3D<CR><LF>

<STX>1H|^&||ABAXIS, INC.^piccolo xpress^2.1.9^0000V0000D|||||D|E 1394-97|20060401160154<CR><ETX>E7<CR><LF>

<STX>2C|1|I|^ ^Received Message NOT Supported|G<CR><ETX>F5<CR><LF>

<STX>3L|1|N<CR><ETX>06<CR><LF>