Table 1: Message Header Record: ASTM Test Tool to Analyzer

ASTM	Field Name	Max Length	Field Contents	Field Description
7.1.1	Record Type ID	1	Н	Header, Upper Case Only
7.1.2	Delimiter Definition	4	I	Field delimiter, vertical bar, (0x7C)
			\	Repeat delimiter, backslash, (0x5C)
			٨	Component delimiter, caret, (0x5E)
			&	Escape delimiter, ampersand, (0x26)

Example of Message Header Record:

H | \(\frac{4}{6} \) \(\text{CR} \)

Table 2: Request Information Record: ASTM Test Tool to Analyzer

ASTM	Field Name	Max Length	Field Contents	Field Description
12.1.1	Record Type ID	1	Q	Query, Upper Case Only
12.1.2	Sequence Number	5	1	A sequence number consistent with sequence number rules. [always 1]
12.1.3	Starting Range ID Number	14	Alphanumeric	Patient ID or Control ID as entered by the Operator on the Analyser
				All Stored Results (Qualified by Time Below) [case insensitive]
			"ALL"	
12.1.7	Beginning Request Results Date and Time	14	YYYYMMDDHHMMS S	Oldest Result that will be transmitted
				If omitted all available will be transmitted
12.1.8	Ending Request Results Date and Time	14	YYYYMMDDHHMMS S	Newest result that will be transmitted
				If omitted all results after beginning time will be transmitted

Examples of Query Record:

Transmit all stored results:

 $Q \mid 1 \mid ALL \langle CR \rangle$

Transmit all stored results after 17 November 2005 at 08:00:00:

Q|1|ALL||||20051117080000<CR>

Transmit all stored results for patient 655321:

Q | 1 | 655321 < CR >

Transmit all stored results for Control Level 1 between 1 January 2005 and 31 December 2005:

Q|1|LEVEL 1||||20050101000000|20051231235959<CR>

Transmit all stored results before 17 November 2005 at 08:00:00:

Q|1|ALL|||||20051117080000<CR>

Table 3: Message Terminator Record: ASTM Test Tool to Analyzer

ASTM	Field Name	Max Length	Field Contents	Field Description
13.1.1	Record Type ID	1	L	Comment, Upper Case Only
13.1.2	Sequence Number	5	1	A sequence number consistent with sequence number rules. [always 1]

Example of Message Terminator Record:

L | 1<CR>